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Nickel-based dental casting alloy intended for metal-ceramic dental restorations without beryllium, cadmium and lead, type 3.

I-BOND 02 is biocompatible non-precious alloy. I-BOND 02 fulfills the recommendations of the standards EN ISO 22674 and EN ISO 9693. The low Vicker's hardness of 180 (HV 10) leads to good and easy milling and grinding behaviour and gives perfect polished surfaces.

Composition	(Mass-%):	Properties
Ni	64,3	Type
Cr	24,2	Dichte
Mo	10	Solidus-, liquidus temperature
Si	1,0	Casting temperature
Co, Fe	< 1 %	Vickers hardness
Nb		WAK
		25 - 500 °C 13,9 x 10 ⁻⁶ K ⁻¹
		20 - 600 °C 14,0 x 10 ⁻⁶ K ⁻¹
		Yield stress
		Rp 0,2 345 MPa
		Modulus of elasticity
		E approx. 202.000 MPa
		Elongation
		A5 26 %

Indications: Used for multiple units of fixed dental prostheses for ceramic firing.

Contraindications: in case of known allergies to alloy content

Intended users: Product is intended to be used by professionals – dental technicians.

Patient selection criteria: Totally or partially edentulous patients.

Recommendations for Use

Wax-up: Do the wax-up as usual. The wall thickness in wax must not be thinner than 0,35 mm. Lead wax sprues indirectly. For wax sprues use round wax wires with Ø 2,5 mm for single crowns and Ø 2,5 – 3,0 mm for bridges. For greater prosthetic constructions with more than 4 pontics use a distribution bar with Ø 4,0 – 4,5 mm.

Investing and Casting: For investing use phosphate bonded investments for crown and bridge work (Intervest K+B Speed or Interfine K+B Speed). Preheat the investment to about 850 to 900 °C. Hold the end temperature of investment for a minimum 30 minutes. Refer to manufacturer's instructions for use for the casting machines.

For I-BOND 02 use an individual ceramic crucible to prevent contamination with other alloys. Clean crucible after every use. When melting by induction heating start casting as soon as the ingots have collapsed and oxide net cracks. For melting by flame heat rotate the reductive zone of the flame around

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