LC INTERTRAY

Version: 04/2026

Date: 8.5.2026

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EN: With the pencil, draw the shape of the tray on a stone model. Use the “Space holder” wax if needed. Insulate the stone model with Isoacryl and let it dry. Put LC Intertray on the model and shape it. Use the cuttings for the holder of the tray. The tray can be polymerised in units that radiate halogen light from a wavelength of 420 nm to 480 nm or UV-light from a wavelength 350 nm to 400 nm. Put the stone model with the tray in the unit and polymerise it for 2 minutes. Remove the tray from the model, turn it upside down and polymerise it again for 2 minutes. Longer polymerisation time improves the strength of the material.

WARNING- Triethylenglykol di-methacrylate, Vinylester Resin

May cause an allergic skin reaction./ Avoid breathing dust, vapours./ Wear protective gloves/protective clothing/eye protection/face protection./ IF ON SKIN: Gently wash with plenty of water and soap./ If skin irritation or rash occurs: Get medical advice/attention./ Take off contaminated clothing and wash it before reuse./ Dispose of contents/container in accordance with local regulations.

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