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EN

Nickel-based dental casting alloy for all-cast crowns and crowns veneered with composite without beryllium, cadmium and lead, type 3.

Non precious nickel-based dental casting alloy for all-cast crowns and crowns veneered with composite. It meets the correspondence to EN ISO 22674.

| Composition | (Mass-%): | Properties |
|-------------|-----------|--|
| Ni | 63,5 | Type 3 |
| Cr | 26 | Density 8,2 g/cm ³ |
| Mo | 8 | Solidus, liquidus temperature 1200°C, 1250°C |
| Si | 1,5 | Casting temperature 1370 °C |
| Fe | 1,4 | Vickers hardness HV 10 190 |
| | | Yield stress Rp 0,2 355 MPa (N/mm ²) |
| | | Modulus of elasticity E cca. 190.000 MPa |
| | | Elongation A5 10 % |

Indications: Used for multiple units fixed dental prostheses with composite covering, not suitable for ceramic firing.

Contraindications: in case of known allergies to alloy content

Intended users: Product is intended to be used by professionals – dental technicians.

Patient selection criteria: Totally or partially edentulous patients.

Recommendations for Use

Wax-up: Do the wax-up as usual. The wall thickness in wax must not be lower than **0,35 mm**. Lead wax sprues indirectly. For wax sprues use round wax wires with **Ø 3,0 mm** for single crowns and **2,5-3,0 mm** for connections between distribution bar and pontics. For greater prosthetic appliances with more than 4 teeth use a distribution bar with **Ø 5,0 mm**.

Investing and Casting: For investing use phosphate bonded investments for crown and bridge work - **Invest K+B Speed** or **Interfine K+B Speed**. Preheat the investment to about **850 to 900 °C**. Hold end temperature of investment for a minimum of 30 minutes. Refer to manufacturer's instructions for use for the casting machines.

For **I-GW** use an individual ceramic crucible to prevent contamination with other alloys. Clean crucible after every use.

When melting by induction heating start casting as soon as the ingots have collapsed and oxide net cracks. For melting by flame heat rotate the reductive zone of the flame around ingots. Start casting as

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soon as the bath begins to vibrate. Allow the cylinder slow air cool down to the ambient temperature and deflask.

Finishing: Sandblast visible surface of the frame with **Al₂O₃ - Interlox 250 µm** and then finish it with suitable grinding and polishing instruments for nickel alloys. At the end use Universal polishing paste for Co-Cr-Mo alloys REF 460 and polish up to high gloss. Clean polished surface with vapour or boiling water.

Soldering and Welding: Soldering before firing of the frame can be carried out with suitable base-metal-solder (**REF 0495 Intersolder**) and high temperature flux (**REF 0496 Interflux**). For welding with laser use suitable base-metal welding wires (**REF 0497 I Weld**). The requirements of EN ISO 9333 must be followed.

Safety warnings: Use suction unit, wear gloves, goggles and protective mask with filter FFP2 when processing the material, to remove and protect yourself against dust.

Reciprocal Actions! In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

Side effects: Consider allergic hypersensitivities to contents of the alloy.

Notice: Any serious incident, that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

Warranty: Technical data based on our own researches and knowledge and can be used as orientation values. The user itself is responsible for appropriate use of the product.

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