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EN

Nickel-based dental casting alloy intended for metal-ceramic dental restorations without beryllium, cadmium and lead, type 3.

I-BOND 02 is biocompatible non-precious alloy. I-BOND 02 fulfills the recommendations of the standards EN ISO 22674 and EN ISO 9693. The low Vickers's hardness of **180 (HV 10)** leads to good and easy milling and grinding behaviour and gives perfect polished surfaces.

Composition	(Mass-%):	Properties
Ni	64,3	Type 3
Cr	24,2	Density 8,2 g/cm ³
Mo	10	Solidus, liquidus temperature 1250 °C, 1340 °C
Si	1,0	Casting temperature 1440 °C
Co, Fe	< 1%	Vickers hardness HV 10 180
Nb		Coefficient of thermal expansion 25 - 500 °C 13,9 x 10 ⁻⁶ K ⁻¹ 20 - 600 °C 14,0 x 10 ⁻⁶ K ⁻¹
		Yield stress Rp 0,2 345 MPa
		Modulus of elasticity E aprox. 202.000 MPa
		Elongation A5 26 %

Indications: Used for multiple units of fixed dental prostheses for ceramic firing.

Contraindications: in case of known allergies to alloy content

Intended users: Product is intended to be used by professionals – dental technicians.

Patient selection criteria: Totally or partially edentulous patients.

Recommendations for Use

Wax-up: Do the wax-up as usual. The wall thickness in wax must not be thinner than **0.35 mm**. Lead wax sprues indirectly. For wax sprues use round wax wires with **Ø 2.5 mm** for single crowns and **Ø 2.5 – 3.0 mm** for bridges. For greater prosthetic constructions with more than 4 pontics use a distribution bar with **Ø 4.0 – 4.5 mm**.

Investing and Casting: For investing use phosphate bonded investments for crown and bridge work (**Intervest K+B Speed** or **Interfine K+B Speed**). Preheat the investment to about **850 to 900 °C**. Hold the end temperature of investment for a minimum **30 minutes**. Refer to manufacturer's instructions for use for the casting machines.

For **I-BOND 02** use an individual ceramic crucible to prevent contamination with other alloys. Clean crucible after every use. When melting by induction heating start casting as soon as the ingots have collapsed and oxide net cracks. For melting by flame heat rotate the reductive zone of the flame around

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ingots. Start casting as soon as the bath begins to vibrate. Allow the cylinder to cool down slowly to the room temperature and then deflask.

Processing: Process the frame in one direction using a hard-metal milling tool. Remove any sharp edges.

Ceramic: Use ceramic, suitable for porcelain to metal restoration. Sand-blast the surface by use of a pencil-blaster.

Blast all of the surface on which the porcelain is to be applied with **Al₂O₃ - Interlox 250 µm** and clean the surface (ultrasonically cleaning, under vapour or boiling water, diethyl acetate). If oxide firing is carried out to check the surface, it need to be fired on **960 °C** at vacuum for **5 - 10 minutes**. After cooling the frame has to show slightly dark greenish colour. Oxide layer should be removed with **Al₂O₃ - Interlox 250 µm**. The frame is cleaned under vapour or boiling water. The opaque is then applied on the surface by a first thin wash firing followed by a second covering opaque. Firing has to be done according to the ceramic manufacturer's instructions. **After each firing step (dentine firing, build-up and glazing) cool down in a normal manner.** Remove ceramic mechanically. Hydrofluoric acid (HF) corrodes the metal.

Finishing: After firing of the metal-ceramic prosthetic work sand-blast the visible surface of the frame with **Al₂O₃ - Interlox 50 µm** and then polish it with suitable grinding and polishing instruments for nickel alloys. At the end use **Universal polishing paste for Co-Cr-Mo alloys REF 460** and polish up to high gloss. Clean polished surface with vapour or boiling water.

Soldering and Welding: Soldering before firing of the frame can be carried out with suitable base-metal-solder and high temperature flux. For welding with laser use suitable base-metal welding wires. The requirements of EN ISO 9333 must be followed.

Safety warnings: Use suction unit, wear gloves, goggles and protective mask with filter FFP2 when processing the material, to remove and protect yourself against dust.

Reciprocal Actions! In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

Side effects: Consider allergic hypersensitivities to contents of the alloy.

Notice: Any serious incident, that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

Warranty: Technical data based on our own researches and knowledge and can be used as orientation values. The user itself is responsible for appropriate use of the product.

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