# CC DISK EASY CoCr

Potrebno upoštevati navodila za uporabo / Consult instructions for use / Bitte beachten Sie die Bedienungsanleitung / Viz návod k použití / Potrebno poštivati uputstva za upotrebu / Consultar las instrucciones de uso / Необходимо учитывать инструкции по применению

Uporabno do / Use by / Gültig bis / Použití / Upotrijebiti do / Usar antes / Использовать до

Kontrolna številka / Batch number / Seriennummer / Seriové číslo / Serijski broj / Numero de serie / Серийный номер

Kataloška številka / Catalogue number / Katalognummer / Katalogové číslo / Kataloški broj / Numero de catalogo / Каталожный номер

Medicinski pripomoček / Medical device / Medizinprodukt / Medicinski proizvod / Zdravotniški pripomoček / Producto sanitaria / Медицинское изделие

Edinstveni identifikator pripomočka / Unique device identifier / Einmalige Produktkennung / Jedinstvena identifikacija proizvoda / Уникальный идентификатор изделия / Jedinečným identifikátorem prostředku / Identificador único del producto

Proizvajalec / Manufacturer / Hersteller / Výrobce / Proizvodač / Fabricante / Производитель

Datum proizvodnje / Date of manufacture / Herstellungsdatum / Datum výroby / Datum proizvodnje / Fecha de fabricacion / Дата изготовления

Navodila za uporabo / Directions for use  
 Gebrauchsanweisung / Instrucciones de uso / Návod k použití  
 Uputa za uporabu / Инструкции по использованию

SR

INTERDENT d.o.o.  
 Zemunska 22, lok 3  
 RS-11070 Novi Beograd  
 T/F: +381/11 217 53 74  
 www.interdent-bg.com

SK

INTERDENT SK s.r.o.  
 Za dráhou 21  
 SK-902 01 Pezinok  
 T: 0903 418 001  
 E: interdent@interdent.sk

CZ

INTERDENT s.r.o.  
 Foerstrova 12, Strašnice  
 CZ-10000 Praha  
 T: +420/274 783 114  
 F: +420/274 820 130  
 E: interdent@interdent.cz

HR

INTERDENT d.o.o.  
 Vinogradski odvojak 2d  
 HR-10431 Sveta Nedelja  
 T: +385/1 3873 644  
 F: +385/1 38736 17  
 E: interdent@interdent.hr

Proizvajalec / Producer / Hersteller  
 Interdent d.o.o. - SI - 3000 CELJE  
 Opekarniška cesta 26 · T: +386 (0)3 425-62-00  
 E: info@interdent.cc · www.interdent.cc

Verzija: 05/2022  
 Datum: 15.12.2022  
 Made in Slovenia

## SL

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## EN

Cobalt-based casted discs, intended for producing prosthetic substitutes in CAD/CAM milling machines for metal-ceramic dental restorations, type 4.

CC DISK EASY CoCr is nickel, beryllium, cadmium and lead free and fulfills the recommendation of the standard EN ISO 22674 and EN ISO 9693 for dental alloys.

The alloy has excellent corrosion resistance and good milling properties. It is easy to drill in and gives perfectly polished surface. Only few oxides are formed during ceramic firing.

Composition (%):	Technical data		
Co 62,5	Type	4	
Cr 27,2	Temperature solidus, liquidus	1380 °C, 1450 °C	
W 8,2	Vickers hardness	HV 10 249	
Si 1,7	Coefficient of thermal expansion	25 - 500 °C 14,4 x 10 <sup>-6</sup> K <sup>-1</sup>	
Mn <1,0		20 - 600 °C 14,6 x 10 <sup>-6</sup> K <sup>-1</sup>	
	Elongation limit	Rp 0,2 380 MPa	
	Tensile strength	Rm 553 MPa	
	Modul of elasticity (GPa)	167 GPa	
	Ductile yield	A5 16,9 %	

**Indications:** Used for multiple units fixed dental prostheses for ceramic firing, also for – for appliances with thin sections that are subject to very high forces, e.g. removable partial dentures, clasps, thin veneered single crowns, full arch fixed dental prostheses or those with small cross-sections, bars, attachments, implant retained superstructures.

**Contraindications:** in case of known allergies to alloy content

**Intended users:** Product is intended to be used by professionals – dental technicians.

**Patient selection criteria:** Totally or partially edentulous patients.

**Ceramic firing:** You can use all types of ceramic which you also use in conventional method on cast crowns and bridges (VITA, Ivoclar, Ducera, etc.).

**Please be careful about cooling procedure. Cobalt-based alloys must be cooled down slowly to 600 °C after each firing step (dentin, correction, glaze).**

1. Remove connectors and smoothen the sharp edges. Sand-blast the surface with 250 microns aluminum oxide on 2-3 bar air pressure (REF 404 Interlox 250 micr).
2. Clean the surface in distilled water with ultrasound, degrease it with ethyl acetate or under vapor.

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5. When "metal-ceramic bond" phase is finished, proceed with opaque layering, to cover the metal shade of the construction.
6. **For firing procedure, consult ceramic manufacturer's instructions for use.**
7. After each firing step (dentine bake, correction bake, glaze firing), cool down slowly to 600 °C.

**Finishing:** After last firing step finish the visible surface of the alloy with polishing rubbers and polishing paste (REF 461 Diamond polishing paste) to high gloss. On the cervical parts of the crowns, remove the oxides with glass pearls or very carefully with 50 microns aluminum oxide (REF 406 Interlox 50 micr.) Be careful not to damage the gingival fit!

**Safety warnings:** Use suction unit, wear gloves, goggles and protective mask with filter FFP2 when processing the material, to remove and protect yourself against dust.

**Reciprocal Actions!** In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

**Side effects:** Consider allergic hypersensitivities to contents of the alloy.

**Notice:** Any serious incident, that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

**Warranty:** Technical data based on our own researches and knowledge and can be used as orientation values. The user itself is responsible for appropriate use of the product.

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