# INTERACRYL ORTHO

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**INTERACRYL ORTHO** (resin for orthodontics)

Interacryl ortho is a self-curing methyl methacrylate for use in orthodontics and is supplied in the polymer/monomer form. The monomer is available in several shades. Interacryl ortho is in accordance with demands of standard EN ISO 20795-2.

**Indications:**

For the removable prosthetics only:

- Orthodontic appliances

**Contraindications:** The possibility of an allergic reaction to PMMA. Unpolymerised or poorly polymerised material don't place in the mouth

**Intended users:** Product is intended to be used by professionals – dental technicians.

**Directions for use:** After thoroughly soaking the plaster model in water, carefully brush on an alginate based separating solution.

**Spray-on technique:** Spray on the first layer of powder, using a jet dispenser bottle. Wet the powder with liquid. Once the powder has soaked up the liquid completely, apply the next layer. Apply alternating layers of powder and liquid to the model. Once the appliance has the desired thickness, we recommend wetting the entire surface with monomer. As soon as the surface loses its sheen, and forms a matt layer, place the appliance and model in the pressure vessel. The water in the vessel should be heated to approx. 45-55 °C. Polymerisation is carried out at a pressure of 2,5 bar, for 20 minutes.

**Dough technique:**

**Mixing ratio at room temperature:**

**Mixing ratio: 14,7 powder / 7ml liquid (2,1 g powder / 1 ml liquid)**

Mix the Interacryl ortho powder and liquid thoroughly for 20 seconds. Allow a doughing time of approx. 95 minutes, at a room temperature of 22 °C, and apply the dough to the model. Assuming a room temperature of approx. 22 °C, the working time will be approx. 6 minutes. The Interacryl ortho appliance must be cured in a pressure vessel at 2.5 bars, for 20 minutes. The water in the vessel should be heated to approx. 45-55 °C.

**Caution:** Do not mix polymer powder with another supplier monomer liquid and vice versa. Place the appliance in water at 50 °C for 1 hour or at room temperature water for 24 hours prior placing in patient mouth.

**Danger:** Highly flammable liquid and vapour. / May cause respiratory irritation. / Causes skin irritation. / May cause an allergic skin reaction. / Causes serious eye irritation. / Keep away from heat/sparks/open flame/hot surfaces. — No smoking. / Keep container tightly closed. / Avoid breathing vapours. / Keep away from heat/sparks/open flames/hot surfaces. — No smoking / Keep container tightly closed. Avoid breathing vapours. / IF ON SKIN: Wash with plenty of soap and water. / Wear protective gloves/protective clothing/eye protection/face protection. / IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing./If eye irritation persists: Get medical advice/attention. / Dispose of contents/ containers in accordance with local regulation.

**Personal protection:** When working with product use protective nitrile gloves, when it is possible to inhale monomer use protective mask with filter A1, when processed polymerized product use protective mask with filter FPP1.

**Side effects:** Possibility of an allergic reaction to MMA and PMMA..

**Notice:** Any serious incident, that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

**Storage:**

- at room temperature
- in the place with good ventilation
- keep away from sunlight
- keep away from children

**Warranty:** Any technical consulting regarding our products, whether verbal, in writing or during demonstrations, is carried out to the best of our knowledge and is to be considered nonbinding. We have no control over use, application and processing, for which the user is responsible.

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