SLO

LOVILEC MAVCA
REF: 620

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EN

PLASTER TRAP
REF: 620


Plaster trap with a closed 4-chamber cleaning system, transparent collector inlaid with foil and fast closing. Transparency tank allows you to visually control the amount of sludge. With a foil input the tank can be quickly and easily cleaned.

1. Technical data:


Height, diameter	320 mm, 40 mm
Drain height	245 mm
Drain diameter	50 mm
Volume	14 l
Supply connection	1½ col, Ø28 mm
Weight	2,1 kg

2. Installation of the plaster trap

1. Remove the U bend from the sink and the associated waste water pipes.
2. Set up the plaster trap in the desired position. The side drain of the plaster trap must be higher than the drain pipe in the wall. The optimal slope should be between 1-2 cm per meter.
3. Mount the inlet hose on the sink. Shorten the inlet hose to the appropriate length so that it can be optimally mounted on the connection on the cover.
4. Make sure that there is no bend upwards on the inlet hose.
5. Mount the drain hose on the side opening and connect it to your sewage pipe. The slope should be 1-2 cm per meter. The supplied clamps can be used so that there is no bend on the sewage hose.

3. Pictorial display of use

Lead the supply tube into the highest chamber.

When cleaning release the screw on the drain and unscrew the locking lid. Skim the standing water and pull out the input.

Remove the pipe-fitting from the

Tie the foil collector under the drain opening and remove it.

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