# INTERACRYL CAST

(SL) Navodila za uporabo / (EN) Directions for use

(DE) Gebrauchsanweisung / (HR) Uputa za korištenje

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**INTERACRYL CAST** (powder and liquid for pouring technique)

Interacryl cast is cadmiumfree autopolymerizable resin for fabricating acrylic dentures. It is supplied as a polymer - powder and monomer - liquid. It is in accordance with demands of EN ISO 20795-1.

**Indications:**

For the removable prosthetics only:

- Partial dentures
- Relines
- Additions / repairs to acrylic dentures
- Pouring acrylic denture bases
- Full dentures
- Implant supported dental prosthesis

**Contraindications:** Conventional method for making partial removable dentures is not possible. The possibility of an allergic reaction to PMMA/MMA. Unpolymerised or poorly polymerised material don't place in the mouth.

**Intended users:** Product is intended to be used by professionals – dental technicians.

**Patient selection criteria:** Totally or partially edentulous patients.

**Directions for use:** Place a silicone or plaster overcast around the teeth on the wax and brush Isoacryl onto the plaster model. Should acrylic teeth have been used, slightly roughen those areas which contact the acrylic base.

**Mixing ratio:**

**colours 3, 5, 15:** 23 g powder (3, 5, 15) / 14 ml liquid for colours 3, 5, 15 (1,64 g powder / 1 ml liquid)

**colour 10:** 23 g powder (10) / 14 ml liquid for color 10 (1,64 g powder / 1 ml liquid)

Pour liquid in mixing bowl and then add powder. Start mixing immediately (30 sec). Ensure that the two materials are mixed thoroughly. In order to take full advantage of the increased working time, one should attempt to attain the desired consistency immediately, without the need to add further powder or liquid. The working time depends on the room temperature (22 °C) and the mixing ratio. The material remains pourable for approx 3.5 minutes after mixing. The doughing time begins approx 4 minutes after mixing and lasts **approx 5 ± 2 minutes (colours 3, 5, 15, 10)**. The material can be kneaded and shaped during this time. The denture can be placed in the pressure vessel when the surface of the acrylic is no longer moist, lose shine on the surface

Water temperature: 55 °C - 60 °C

Operating pressure: 2.5 bar

Polimerization time: 35 - 50 min (depends on thickness of the denture)

**Caution:** Do not mix polymer powder with another supplier monomer liquid and vice versa. In order to ensure that the denture cures properly, it must not be allowed to contact other acrylics, e.g. remains of powder and liquid in measuring beakers, mixing vessels and on spatulas. Place the appliance in water at 50 °C for 1 hour or at room temperature water for 24 hours prior placing in patient mouth.

**Danger:** Highly flammable liquid and vapour. / May cause respiratory irritation. / Causes skin irritation. / May cause an allergic skin reaction. / Causes serious eye irritation. / Keep away from heat/sparks/open flame/hot surfaces. - No smoking. / Keep container tightly closed. / Avoid breathing vapours. / Keep away from heat/sparks/open flames/hot surfaces. - No smoking / Keep container tightly closed. Avoid breathing vapours. / IF ON SKIN: Wash with plenty of soap and water. / Wear protective gloves/protective clothing/eye protection/face protection. / IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing./ If eye irritation persists: Get medical advice/attention. / Dispose of contents/ containers in accordance with local regulation.

**Personal protection:** When working with product use protective nitrile gloves, when it is possible to inhale monomer use protective mask with filter A1, when processed polymerized product use protective mask with filter FFP1.

**Side effects:** Possibility of an allergic reaction to MMA and PMMA.

**Storage:**

- at room temperature
- in the place with good ventilation
- keep away from sunlight
- keep away from children

**Notice:** Any serious incident, that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

**Warranty:** Any technical consulting regarding our products, wheather verbal, in writing or during demonstrations, is carried out to the best of our knowledge and is to be considered nonbinding. We have no control over use, application and processing, for which the user is responsible.

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