INTERACRYL COLD

- (SL) Navodila za uporabo / (EN) Directions for use
(DE) Gebrauchsanweisung / (HR) Uputa za korištenje
(RU) Инструкции по использованию / (MK) Препорака за употреба
(CS) Návod na zpracování / (ES) Instrucciones de uso
(KK) Пайдалану бойынша нұсқаулық / (GR) Διαδικασία επεξεργασίας
(RO) Instrucțiuni de utilizare / (ALB) Pergatitja / (FR) Instructions de procédure

CS

INTERDENT s.r.o.
Foerstrova 12, Strašnice
CZ-10000 Praga
T: +420/274 783 114
F: +420/274 820 130
E:interdent@interdent.cz

HR

INTERDENT d.o.o
Vinogradski odvojak 2d
HR-10431 Sveta Nedelja
T: +385/1 3873 644
F: +385/1 38736 17
E: interdent@interdent.hr

SK

INTERDENT SK S.R.O.
Za dráhou 21
SK-902 01 Pezinok
T: 0903 418 001
E: interdent@interdent.sk

SR

INTERDENT d.o.o.
Egipatska 22, lok 3
RS-11070 Novi Beograd
T/F: +381/11 217 53 74
www.interdent-bg.com

Proizvajalec/Producer/Hersteller
Interdent d.o.o. · SI - 3000 CELJE · Opekarniška cesta 26
T: +386 (0)3 425-62-00 · E: info@interdent.cc
www.interdent.cc

Type: 2
Class:1
CE 0197

Verzija: 14/2025
Datum: 4.2.2025

Made in Slovenia

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INTERACRYL COLD

Interacryl cold is a cadmium-free, methyl methacrylate repair acrylic, supplied in the polymer/monomer form. It contains a catalyst which ensures fast, uniform curing. The colour of the acrylics corresponds with the colour of other acrylics of the production in INTERDENT. Interacryl cold is in accordance with demands of standard EN ISO 20795-1.

Indications:

- Repairing full and partial dentures
- Relining full and partial dentures
- Extending dentures

Contraindications: Not for full total or partial dentures. The possibility of an allergic reaction to PMMA/MMA. Unpolymerised or poorly polymerised material don't place in the mouth.

Intended users: Product is intended to be used by professionals – dental technicians.

Patient selection criteria: Totally or partially edentulous patients.

Directions for use: After carefully removing all wax from plaster model, brush Isoacryl mould seal onto it. The edges of the repair must be roughened and moistened with monomer, prior to applying Interacryl cold.

Mixing ratio at room temperature:

colours 3, 5, 15: 12 g powder (3, 5, 15) / 7 ml liquid for colours 3, 5, 15 (1,7 g powder / 1 ml liquid)

color 10, transparent: 13 g powder (10, transparent) / 10 ml liquid for colours 10 and transparent (1,3 g powder / 1 ml liquid)

Pour liquid in mixing bowl and then add powder. Start mixing immediately. Mix the powder and liquid thoroughly for 20 seconds. Any air bubbles should be eliminated by knocking on the mixing vessel. Allow the mixture a few minutes doughing time, until the desired consistency has been reached. Assuming a room temperature of 22 °C, the working time will be aprox 5 minutes. We recomend curing the repair in a pressure vessel at 2.5 ba for 15 minutes. The water should be heated to 55 - 60 °C.

Doughing time: not defined

Working time (at 22 °C): 5 minutes

Curing time (55 - 60 °C): 15 minutes (under pressure)

After removing the plaster, we process polimerized object with appropriate instruments and polish with polishing preparations (pumice, interlesk, polishing rubbers and brushes).

Caution: Do not mix polymer powder with another supplier monomer liquid and vice versa. Place the denture in water at 50 °C for 1 hour or at room temperature water for 24 hours prior placing in patient mouth.

Danger: Highly flammable liquid and vapour. / May cause respiratory irritation. / Causes skin irritation. / May cause an allergic skin reaction. / Causes serious eye irritation. / Keep away from heat/sparks/open flame/hot surfaces. - No smoking. / Keep container tightly closed. / Avoid breathing vapours. / Keep away from heat/sparks/open flames/hot surfaces. - No smoking / Keep container tightly closed. Avoid breathing vapours. / IF ON SKIN: Wash with plenty of soap and water. / Wear protective gloves/protective clothing/eye protection/face protection. / IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing./ If eye irritation persists: Get medical advice/attention. / Dispose of contents/containers in accordance with local regulation.

Personal protection: When working with product use protective nitrile gloves, when it is possible to inhale monomer use protective mask with filter A1, when processed polimerized product use protective mask with filter FFP1.

Side effects: Possibility of an allergic reaction to MMA and PMMA.

Storage:

- at room temperature
- in the place with good ventilation
- keep away from sunlight
- keep away from children

Notice: Any serious incident, that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

Warranty: Any technical consulting regarding our products, wheather verbal, in writing or during demonstrations, is carried out to the best of our knowledge and is to be considered nonbinding. We have no control over use, application and processing, for which the user is responsible.

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