Isoprint

Water-based, alcohol-free isolator for 3D printed models

Interdent d.o.o. · SI - 3000 CELJE · Opekarniška cesta 26
T: +386 (0)3 425-62-00 · E: info@interdent.cc · www.interdent.cc

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EN

Usage: ISOPRINT is a water-based, alcohol-free isolator intended for the isolation of SLA/DLP 3D-printed models when applying cold polymerising acrylic materials. It provides good insulation and is easily removed with warm water.

Application: For good insulation, apply Isoprint in two layers. Squeeze the necessary amount of Isoprint from the bottle and spread it in a thin layer using a brush (e.g. REF A01).

Dry the applied layer:

1. using a hairdryer for approx. 1 min. or
2. under a polymerisation lamp at T 45-50°C for 3-4 min. or
3. air dry for 20 min.

Before applying the second layer, the first must be completely dry so that it is not sticky to the touch.

Before insulating, block any undercuts on the model with wax (e.g. REF 186-30 Block-Out Wax). After the insulation, continue with the application of the acrylic mass according to the manufacturer's instructions.

After polymerisation, lift the polymerised acrylic product from the printed model with the help of a carving knife. Remove the remains of the insulating agent from the model and acrylic using warm water with the help of a brush or steam with a steam cleaner.

After use, close the bottle (rotate the applicator clockwise) to prevent evaporation.

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