Potrebno upoštevati navodila za uporabo / Consult instructions for use / Bitte beachten Sie die Bedienungsanleitung / Viz návod k použití / Potrebno poštivati uputstva za upotrebu / Consultar las instrucciones de uso / Необходимо учитывать инструкции по применению

Uporabno do / Use by / Gültig bis / Použití / Upotrijebiti do / Usar antes / Использовать до

Kontrolna številka / Batch number / Seriennummer / Seriové číslo / Serijski broj / Numero de serie / Серийный номер

Kataloška številka / Catalogue number / Katalognummer / Katalogové číslo / Kataloški broj / Numero de catalogo / Каталожный номер

Medicinski pripomoček / Medical device / Medizinprodukt / Medicinski proizvod / Zdravotniški pripomoček / Producto sanitaria / Медицинское изделие

Edinstveni identifikator pripomočka / Unique device identifier / Einmalige Produktkennung / Jedinствена identifikacija proizvoda / Уникальный идентификатор изделия / Jedinečným identifikátorem prostředku / Identificador único del producto

Proizvajalec / Manufacturer / Hersteller / Výrobce / Proizvodač / Fabricante / Производитель

Datum proizvodnje / Date of manufacture / Herstellungsdatum / Datum výroby / Datum proizvodnje / Fecha de fabricacion / Дата изготовления

**SL** Navodila za uporabo / **EN** Directions for use  
**DE** Gebrauchsanweisung / **ES** Instrucciones de uso / **CS** Návod k použití  
**HR** Uputa za uporabu / **RU** Инструкции по использованию

**SR**

INTERDENT d.o.o.  
Egipatska 22, lok 3  
RS-11070 Novi Beograd  
T/F: +381/11 217 53 74  
www.interdent-bg.com

**SK**

INTERDENT SK s.r.o.  
Za dráhou 21  
SK-902 01 Pezínok  
T: 0903 418 001  
E: interdent@interdent.sk

**CS**

INTERDENT s.r.o.  
Foerstrova 12, Strašnice  
CZ-10000 Praga  
T: +420/274 783 114  
F: +420/274 820 130  
E:interdent@interdent.cz

**HR**

INTERDENT d.o.o.  
Vinogradski odvojak 2d  
HR-10431 Sveta Nedelja  
T: +385/1 3873 644  
F: +385/1 38736 17  
E: interdent@interdent.hr

Proizvajalec / Producer / Hersteller  
Interdent d.o.o. · SI - 3000 CELJE  
Opekarniška cesta 26 · T: +386 (0)3 425-62-00  
E: info@interdent.cc · www.interdent.cc

**CE** 0197

Verzija: 12/2025

Datum: 29.01.2025

Made in Slovenia

**SL**

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**EN**

Cobalt-based casted discs, intended for making prosthetic substitutes in CAD/CAM milling machines for metal-ceramic dental restorations, type 4.

CC DISK NF CoCr is nickel, beryllium, cadmium and lead free and fulfills the recommendation of the standard EN ISO 22674 and EN ISO 9693 for dental alloys. Vickers hardness 281 (HV10) value leads to good milling and drilling behavior and gives perfect polished surface. Only few oxides are formed during ceramic firing.

Composition	(Mass %):	Properties		
Co	63	Type		4
Cr	24	Density		8,3 g/cm <sup>3</sup>
W	8	Vickers hardness	HV 10	281
Mo	3	Coefficient of thermal expansion	25 - 500 °C	13,9 x 10 <sup>-6</sup> K <sup>-1</sup>
Si	1,0		20 - 600 °C	14,0 x 10 <sup>-6</sup> K <sup>-1</sup>
Nb	<1%	Yield stress	Rp 0,2	478 MPa
		Tensile strength	Rm	581,5 MPa
		Modulus of elasticity	E	approx. 210.000 MPa
		Elongation	A5	10,4 %

**Indications:** Used for multiple units fixed dental prostheses for ceramic firing, also for – for appliances with thin sections that are subject to very high forces, e.g. removable partial dentures, clasps, thin veneered single crowns, full arch fixed dental prostheses or those with small cross-sections, bars, attachments, implant retained superstructures.

**Contraindications:** in case of known allergies to alloy content

**Intended users:** Product is intended to be used by professionals – dental technicians.

**Patient selection criteria:** Totally or partially edentulous patients.

**Ceramic firing:** You can use all types of ceramic which you also use in conventional method on cast crowns and bridges (VITA, Ivoclar, Ducera, etc.).

**Ceramic firing:** You can use all types of ceramic which you also use in conventional method on cast crowns and bridges (VITA, Ivoclar, Ducera, etc.).

**Please be careful about cooling procedure, Co- based alloys must be cooled down slowly after each firing step.**

- Remove connectors and smooth sharp edges. Sand blast the surface with 150 or 250 microns

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- aluminum oxide (Interlox 250 micr or Interlox 150 micr).
- Clean the surface in distilled water with ultrasound, degrease with ethyl acetate or under vapor.
- Carry out oxide firing (only to control the surface) 5 min at 960 °C. After oxide firing removing of oxides is required. Use aluminum oxide 250 or 150 micr and again degrease the surface.
- First layer, opaque, applied thinly and evenly over the construction, then follow direction for use of ceramic producer.
- Second layer of the opaque apply the way that all the surface of construction is covered and that the metal shade is not visible over the ceramic.
- After each firing carry out slowly cooling down.

**Finishing:** After last firing step finish visible surface of the alloy with polishing rubbers and polishing paste for Co-Cr alloys (REF 0460 Universal polishing paste for Co-Cr) to high gloss. Remove oxides with glass pearls or very carefully with aluminum oxide 50 microns on the cervical parts of the crowns (that gingival fit is not damaged).

**Safety warnings:** Use suction unit, wear gloves, goggles and protective mask with filter FFP2 when processing the material, to remove and protect yourself against dust.

**Reciprocal Actions!** In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

**Side effects:** Consider allergic hypersensitivities to contents of the alloy.

**Note:** Any serious incident, that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

**Warranty:** Technical data based on our own researches and knowledge and can be used as orientation values. The user itself is responsible for appropriate use of the product.

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