KORENINSKI ZATIČI ROOT CANAL POSTS

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GB ROOT CANAL POSTS

Root posts are made out of stainless steel. They are used in endodontically treated teeth that has not have enough core material. Treatment with root post can be successful only if there is 1,5-2 mm of dentine ferrule that provide support of root post. If dentine ferrule is not provided than root post can be used as temporary solution.

Root posts has to be disinfected, cleaned, dried and sterilized before use. Sterilize it in autoclave with steam at temperature 134 °C, pressure 2,1 bar for 5 minutes.

- Treat the root of the tooth endodontically.
- Remove as much filling of the root as possible.
- Prepare the root canal with the drill.
- Make sure about the length of the canal by measurement or X-ray scann.
- Drill desired length in the root. Recommended velocity: approx. 500-1000 rpm.
- To test the preparation of the canal, place the root post in it.
- Clean and disinfect the root canal.
- Fill in the canal with cement.
- Wet the area of the post that will be placed in the root.

- Place the post slowly into the root canal.
- Shape the upper part of the root with composite or temporary bridge and crown material.
- Take an impression for making ceramic or acrylic crown.

Warning! In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur which lead to not predictable solubility of metal.

Side effects: Allergic reaction to the composition of the root post is possible. Dentist chooses appropriate material.

Warranty! Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our products are subject to further development. Therefore alterations in construction and composition are reserved.

Element	Composition mass %	Characteristic:
Fe	63	Elongation limit Rp 0.2 1100 MPa
Cr	17,4	Tensile strength Rm 1600 MPa
Ni	14,4	Ductile yield A5 5 %
Mo	2,71	
Mn	1,78	
Si	≤ 1	
C, P, N	< 1	

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