# I-MG

<b>CS</b>	<b>HR</b>	<b>SR</b>	<b>SK</b>
INTERDENT s.r.o. <div>Foerstrova 12, Strašnice CZ-10000 Praha T: +420/274 783 114 F: +420/274 820 130 E:interdent@interdent.cz</div>	INTERDENT d.o.o. <div>Vinogradski odvojak 2d HR-10431 Sveta Nedelja T: +385/1 3873 644 F: +385/1 3873617 E: interdent@interdent.hr</div>	INTERDENT d.o.o. <div>Egipatska 22, Iok 3 RS-11070 Novi Beograd T / F: +381/11 217 53 74 www.interdent-bg.com</div>	INTERDENT SK s.r.o. <div>Za dráhou 21 SK-902 01 Pezinok T: 0903 418 001 E: interdent@interdent.sk</div>

<b><span><span>Interdent</span></span></b> <sup>®</sup> 	<b>CE 0197</b>
Interdent d.o.o. · SI · 3000 CELJE <div>Opekarniška cesta 26 · T: +386 (0)3 425-62-00 <div>E: info@interdent.cc · <b>www.interdent.cc</b></div></div>	Version: 18/2025 <div>Date: 25.11.2025 <div>Made in Slovenia</div></div>

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<b>EN</b>
-----------

Directions for use

**Cobalt-based dental casting alloy for partial dentures without beryllium, nickel, cadmium and lead, type 5.**
**I-MG** is a biocompatible alloy. It has outstanding resistance to corrosion and a high strength. **I-MG** can be polished and laser-welded. It meets the requirements of EN ISO 22674.

Composition	m <span> </span> %	Technical data	
<b>Co</b>	62,5	<b>Type</b>	5
<b>Cr</b>	29,5	<b>Density</b>	8,2 g/cm <sup>3</sup>
<b>Mo</b>	5,5	<b>Solidus, liquidus temperature</b>	1295 <span> </span> °C, 1345 <span> </span> °C
<b>Si</b>	1,4	<b>Casting temperature</b>	1450 <span> </span> °C
<b>Mn, C,N</b>	<1,0	<b>Vickers hardness</b>	HV 10 365
		<b>Tensile strength</b>	R <sub>m</sub> 914 MPa
		<b>Yield stress</b>	R <sub>p0,2</sub> 640 MPa
		<b>Modulus of elasticity</b>	E ≈220 GPa
		<b>Elongation after fracture</b>	A <sub>5mm</sub> 7,5%

**Indications:** Used for frames in removable prosthetic dentures, for appliances in which parts require the combination of high stiffness and stress resistance, e.g. thin removable partial dentures, parts with thin cross-sections, clasps, attachments, bars.
**Contraindications:** In case of known allergies to alloy content.
**Intended users:** The product is intended to be used by professionals – dental technicians.
**Patient selection criteria:** Tolly or partially edentulous patients.
**Recommendations for use**
**Wax-up:** Sprues should have a round profile with Ø 3.5–4 mm, avoid rectangular placing of sprues and direct slap of alloy.
**Investing:** Use phosphate-bonded partial-denture investment material REF 933- Modelcast or REF 933S- Modelcast S. Preheating temperature of investment is 950–1000°C and depends on the moulding and structure of sprues. The end temperature at total for plates is 1050 °C. The holding time of the end temperature is 45–60 minutes, and depends on the size and the number of the casting flask. Please follow the recommendations for use of the casting machine manufacturer for casting.

**Casting:** For **I-MG**, use a new ceramic crucible to prevent contamination with other alloys. Clean the crucible after every use. When melting by induction heating, start casting as soon as the ingots have collapsed and the oxide net cracks. For melting by flame heat, rotate the reductive zone of the flame around asgones. Start casting as soon as the bath begins to vibrate. Do not overheat the alloy. After cooling down to ambient temperature, deflask the cast and sand blast it under 2–4 bar of pressure with Al<sub>2</sub>O<sub>3</sub> 110–250 µm (Interalox). Be careful with clasp sides and retention elements.

**Polishing:** For finishing, use carbide, ceramically bonded stones or carbide bur tools. Use a polishing liquid for electrolyte polishing (REF 490- Elektrolyt MG) in a dental laboratory. Clasps and fitting parts should be protected before polishing with a special varnish that prevents uncontrolled material removal. After treating with stones, milling tools and measuring the fit, the cast should be polished with a rubber polisher, and at the end with Universal polishing paste for Cr-Co-Mo alloys REF 460 to high gloss.

**Soldering and welding:** Soldering with a Co-based lot (REF 495- Intersolder) and a high-temperature flux (REF 496- Interflux). Laser welding with suitable base-metal welding wires. The requirements of EN ISO 9333 must be followed.

**Safety warnings:** Use a sun protection unit, wear gloves, googles and a protective mask with an FFP2 filter when processing the material to remove and protect yourself against dust.

**Reciprocal actions:** In case of occlusal or approximal contact of different alloys, electrochemical reactions may occur in very rare instances.

**Side effects:** Consider allergic hypersensitivities to the contents of the alloy.

**Notice:** Any serious incident that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

**Warranty:** Whether given verbally, in writing, or by practical instructions, our recommendations for use are based upon our own experience and trials, and can only be considered as reference values. Our products are subject to further development. Therefore, alterations in the structure and composition are reserved.

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