WAX CENTER

REF 0604

 Navodila za uporabo /  Directions for use
 Upute za upotrebu /  Gebrauchsinformationen
 Инструкции по применению

CZ
INTERDENT s.r.o.
Foerstrova 12, Strašnice
CZ-10000 Praha
T: +420/274 783 114
F: +420/274 820 130
E:interdent@interdent.cz

HR
INTERDENT d.o.o.
Vinogradski odvojak 2d
HR-10431 Sveta Nedelja
T: +385/1 3873 644
F: +385/1 38736 17
E: interdent@interdent.hr

SRB
INTERDENT d.o.o.
Zemunska 22, lok 3
RS-11070 Novi Beograd
T/F: +381/11 217 53 74

MK
INTERDENT d.o.o.
Bojmija 1-12M
MK-1000 Skopje
T: +389/224 628 55
F: +389/224 634 12
E: intsk@t-home.mk

SK
INTERDENT SK, s.r.o.
Lišcie údolie 57
SK-84231 Bratislava
T: 00421 2 5440554,
F: 00421 2.5440555,
E: interdent@interdent.sk

INTERDENT®

Proizvajalec/Producer/Hersteller

SI - 3000 CELJE · Opekarniška cesta 26
T: +386 (0)3 425-62-00 · F: +386 (0)3 425-62-10
E: info@interdent.cc · www.interdent.cc

Verzija: 02/2015

Made in Slovenia

WAX CENTER

REF 0604

SI

Kombiniran aparat za modeliranje in potapljanje

NAVODILO ZA UPORABO

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REF 0604

GB

Combined device for dipping and waxing up

DIRECTION FOR USE

1. Usage area

Wax center REF 0604 is combined electrical device for waxing up and dipping method in dental laboratory.

2. Advantages

Suitable for laboratories with two working places

3. Technical data

Connection:	230 V, 50/60 Hz
Fuse:	0.16 A
Power taken away at max. temperature:	20 W
Heating pin:	voltage 6 V/5 W, cable 1.2 m
Temperature range of heated pin:	30 - 230°C without regulation
Temperature range of the pot:	60 - 105°C without regulation
Function signalization:	lightened switcher
Outside measurements:	width = 160 mm, depth = 80 mm, height = 100 mm
Weight:	1,4 kg
Cable with protective connector:	1,5 m

4. Description

Wax Center is a compact, efficient device for waxing up with heated tip and for making a cap with a dipping method. Mechanically adjustable operating temperature enables waxing up tension free and accurate results at dipping method. Temperature of heated pin comprises a range of 30 - 230°C, and the dipping pot the range of 60 - 105°C. Each subsequent reservoir has around 5 - 7°C lower heating temperature in the pot as the previous one. You can not regulate the temperature of each pot separately. Cork sleeve insulates the handpiece well which ensures safe modeling.

4.1. Before first use

When you take new device from the box make sure it is not damaged from the transport. If you notice any damage return the device immediately. Put the device on desired place and connected it with the plug on heating container. Make sure that the net voltage is in accordance with data on the board of device. If it is appropriate than connect the device on the net.

4.2. Usage

4.2.1. Direction

It is very easy to change wax-up instruments. You can do it by hand when the thermo knife is cold and turned off. During work you can change the instruments by special tongs resistant to temperature 250°C. Be careful on burns! Turn on the device with switch. It is ready when there is a green light on the switch. The desired working temperature is placed with swivel knob. Heating point is than prepared for use. Fill the pot with suitable wax. Switch the device on. Green light on the switcher shows that the device is ready. Regulate the desired temperature with swivel knob. Be careful! Touching the edge of the pot can cause burns.

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4.2.2. Mistakes removing

If the pot and waxing up instrument do not heat, if electronic does not work and if even when changing the pin, the device does not give desired results, you must contact authorized service that will find and remove mistake.

4.3. Maintenance

It is very important to maintain clean waxing up instrument and the pot. Repairing of the device must be done by authorized service. Interdent is not responsible for unauthorized intervention in device. Guarantee for device is one year and does not valid in case of using the device in other purposes that are described in this directions for use and in case of intervention in device by unauthorized service.

5. Safety precaution

Danger of burns at direct contact of hot pin and hot area of the pot. When heating pin is not used, it must be lay down on holder to prevent overheating and possible fire. The intervention in device can be made by qualified person. Do not use device on people.

6. Content

REF 0604 Combined device for dipping and waxing up	1 pcs
REF 0606 Waxing up instrument, knife	1 pcs
REF 0607 Waxing up instrument awl	1 pcs
REF 0608 Waxing up instrument, pin, thin 0,5 mm	1 pcs
REF 0609 Waxing up instrument, pin middle 1,0 mm	1 pcs
REF 0617 Cork holder	2 pcs

7. Accessories

REF 0610 Waxing up instrument, knife small,
REF 0611 Waxing up instrument, knife small, transversal
REF 0612 Waxing up instrument, with hollow
REF 0613 Waxing up instrument, opener left, 0,5 mm
REF 0614 Waxing up instrument, opener left, 1,0 mm
REF 0615 Waxing up instrument, pin, extra thin hart gold plated 0,3 mm
REF 0618 Cork holder

8. Warranty

Our recommendations and instructions are based on our own experiences. Using of device in any other purposes as described above is made on own responsibility. Quality of device is tested and assured.

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Назив izdelka: Kombiniran aparat za modeliranje in potopno metodo

Name of the product: *combined device for dipping and waxing up*

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The product will work flawless during the warranty period if being administered in accordance with its purpose and the instructions provided

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Warranty period expires after 12 months from purchase day. At your request we will repair failure or defects on the product, if you will report it within the warranty period. Failure will be resolved free of charge within 45 days from the date of failure report. The product, which will not be repaired within this deadline, will be, at your request, replaced with a new one. We will extend the warranty period for the time the product was being repaired.

- Veljavnost garancije dokazujete s potrjenim garancijskim listom in original računom.

Validity of the warranty should be proofed by certified warranty card and original invoice

Garancija preneha veljati:

- Če v aparat posega nepooblaščen a oseba / If an unauthorized person interferes into the device
- Če niso upoštevana priložena navodila za uporabo / If provided instructions for use are not being considered
- Če je izdelek mehansko poškodovan / If the device is being mechanically injured
- Če so vgrajeni neoriginalni rezervni deli / If the unoriginal spare parts are being build in

SI - 3000 CELJE · Opekarniška cesta 26 · T: +386 (0)3 425-62-00
F: +386 (0)3 425-62-10 · E:info@interdent.cc · www.interdent.cc