Potrebno upoštevati navodila za uporabo / Consult instructions for use / Bitte beachten Sie die Bedienungsanleitung / Potrebno poštivati uputstva za upotrebu / Vodete se по uputstvoto за употреба / Необходимо учитывать инструкции по применению / Viz návod k použití

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GB

Nickel-based dental casting alloy for all-cast crowns and crowns veneered with composite without beryllium, cadmium and lead, type 3.

Non precious nickel-based dental casting alloy for all-cast crowns and crowns veneered with composite. Contains no Beryllium, Cadmium and lead. Not for metal to ceramic restoration. It meets the correspondence to EN ISO 22674.

Composition	(Mass-%):	Properties		
Ni	63,5	Type		3
Cr	26	Density		8,2 g/cm <sup>3</sup>
Mo	8	Solidus, liquidus temperature		1200 °C, 1250 °C
Si	1,5	Casting temperature		1370 °C
Fe	1,4	Vickers hardness	HV 10	190
		Yield stress	Rp 0,2	355 MPa (N/mm <sup>2</sup> )
		Modulus of elasticity	E	cca. 190.000 MPa
		Elongation	A5	10 %

#### Recommendations for Use

##### Wax-up

Do the wax-up as usual. The wall thickness in wax must not be lower than **0,35 mm**. Lead wax sprues indirectly. For wax sprues use round wax wires with **Ø 3,0 mm** for single crowns and **2,5-3,0 mm** for connections between distribution bar and pontics. For greater prosthetic appliances with more than 4 teeth use a distribution bar with **Ø 5,0 mm**.

##### Investing and Casting

For investing use phosphate bonded investments for crown and bridge work - **Intervest K+B Speed or Interfine K+B Speed**. Preheat the investment to about **850 to 900 °C**. Hold end temperature of investment for a minimum of 30 minutes. Refer to manufacturer's instructions for use for the casting machines.

For I-GW use an individual ceramic crucible to prevent contamination with other alloys. Clean crucible after every use. When melting by induction heating start casting as soon as the ingots have collapsed and oxide net cracks. For melting by flame heat rotate the reductive zone of the flame around ingots. Start casting as soon as the bath begins to vibrate. Allow the cylinder slow air cool down to the ambient temperature and deflask.

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#### Finishing

Sandblast visible surface of the frame with **Al<sub>2</sub>O<sub>3</sub> - Interlox 250 µm** and then finish it with suitable grinding and polishing instruments for nickel alloys. At the end use Universal polishing paste for Co-Cr-Mo alloys REF 460 and polish up to high gloss. Clean polished surface with vapour or boiling water.

#### Soldering and Welding

Soldering before firing of the frame can be carried out with suitable base-metal-solder (REF 0495 Intersolder) and high temperature flux (REF 0496 Interflux). For welding with laser use suitable base-metal welding wires (REF 0497 I Weld). The requirements of EN ISO 9333 must be followed.

#### Final Safety Notes!

Metal dusts are harmful. Use dust extractor, wear gloves, googles and protective mask with filter FFP2 during handling with alloy.

#### Side effects:

Consider allergic hypersensitivities to contents of the alloy.

#### Reciprocal Actions!

In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

#### Warranty!

Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our product are subject to further development. Therefore alterations in construction and composition are reserved.

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