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**Cobalt-based dental casting alloy for partial dentures without beryllium, nickel, cadmium and lead, typ 5.**

**I-MG** is biocompatible alloy. **I-MG** can be polished and laser welded. It meets the requirements of standard EN ISO 22674.

Composition	(Mass-%):	Properties	
Co	62,5	Type	5
Cr	29,5	Density	8,2 g/cm³
Mo	5,5	Solidus, liquidus temperature	1295 °C, 1345 °C
Si	1,4	Casting temperature	1450 °C
Mn, C, N	<1%	Vickers hardness	HV 10 365
		Tensile strength	Rm 914 MPa (N/mm²)
		Yield stress	Rp 0,2 640 MPa (N/mm²)
		Modulus of elasticity	E cca. 220.000 MPa
		Elongation	A5 7,5 <span> </span> %

**Indications:** Used for frames in removable prosthetic dentures, for appliances in which parts require the combination of high stiffness and proof stress, e.g. thin removable partial dentures, parts with thin cross-sections, clasps, attachments, bars.

**Contraindications:** in case of known allergies to alloy content

**Patient selection criteria:** Totally or partially edentulous patients.

**Intended users:** Product is intended to be used by professionals – dental technicians.

**Recommendation for use**

**Wax up:** Sprues should be with round profile with **Ø 3.5 – 4 mm**, avoid rectangular placing of sprues and direct slap of alloy.

**Investing:** Use phosphate-bonded partial-denture investment material Modelcast or Modelcast S. Preheating temperature of investment is **950 – 1000 °C** and depend on moulding and construction of sprues. At total plates is the end temperature **1050 °C**. Holding time of end temperature is

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**45 - 60 minutes**, and depends on the size of the casting flask and filler of the furnace. We ask you to follow recommendation for use of casting machine manufacturer for casting.

**Casting:** For **I-MG** use an individual ceramic crucible to prevent contamination with other alloys. Clean crucible after every use. When melting by induction heating start casting as soon as the ingots have collapsed and oxide net cracks. For melting by flame heat rotate the reductive zone of the flame around ingots. Start casting as soon as the bath begins to vibrate. Do not overheat the alloy. After cooling down to ambient temperature, deflask the cast and sand-blasted it under 4 bar pressure with **Al<sub>2</sub>O<sub>3</sub> - Interalex 250 µm**. Be careful with clasp sides and stress breakers.

**Polishing:** For finishing use carbide, ceramically bonded stones or hart metal milling tools.

Use polishing liquid for electrolyte polishing in dental laboratory. Clasps and fitting parts should be protected before polishing with special varnish that prevent uncontrolled material removing.

After treating with stones, milling tools and measuring the fit, the cast should be polished with rubber polisher, and at the end with **Universal polishing paste for Cr-Co-Mo alloys REF 460** up to high gloss.

**Soldering and welding:** Soldering with Co based lot and high temperature flux. Laser welding with suitable base-metal welding wires. The requirements of EN ISO 9333 must be followed.

**Safety warnings:** Use suction unit, wear gloves, googles and protective mask with filter FFP2 when processing the material, to remove and protect yourself against dust.

**Reciprocal Actions!** In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

**Side effects:** Consider allergic hypersensitivities to contents of the alloy.

**Notice:** Any serious incident, that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

**Warranty!** Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our product are subject to further development. Therefore alterations in construction and composition are reserved.

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