INTERACRYL COLD

(SL) Navodila za uporabo / (EN) Directions for use

(DE) Gebrauchsanweisung / (HR) Uputa za korištenje

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INTERACRYL COLD

Interacryl cold is a cadmium-free, methyl methacrylate repair acrylic, supplied in the polymer/monomer form. It contains a catalyst which ensures fast, uniform curing. The colour of the acrylics corresponds with the colour of other acrylics of the production in INTERDENT. Interacryl cold is in accordance with demands of standard EN ISO 20795-1.

Indications:

- Repairing full and partial dentures
- Relining full and partial dentures
- Extending dentures

Contraindications: Not for full total or partial dentures. The possibility of an allergic reaction to PMMA. Unpolymerised or poorly polymerised material don't place in the mouth.

Intended users: Product is intended to be used by professionals – dental technicians.

Patient selection criteria: Totally or partially edentulous patients.

Directions for use: After carefully removing all wax from plaster model, brush Izoacryl mould seal onto it. The edges of the repair must be roughened and moistened with monomer, prior to applying Interacryl cold.

Mixing ratio at room temperature:

colours 3, 5, 15, 16: 12 g powder / 7 ml liquid (1,7 g powder / 1 ml liquid)
color 10: 13 g powder / 10 ml liquid (1,3 g powder / 1 ml liquid)

Pour liquid in mixing bowl and then add powder. Start mixing immediately. Mix the powder and liquid thoroughly for 20 seconds. Any air bubbles should be eliminated by knocking on the mixing vessel. Allow the mixture a few minutes doughing time, until the desired consistency has been reached. Assuming a room temperature of 22 °C, the working time will be aprox 5 minutes. We recommend curing the repair in a pressure vessel at 2.5 ba for 10 minutes. The water should be heated to 45-50 °C.

Doughing time: not defined

Working time (at 22 °C): 5 minutes

Curing time: 10 minutes (under pressure)

After removing the plaster, we process polymerized object with appropriate instruments and polish with polishing preparations (pumice, interlesk, polishing rubbers and brushes).

Caution: Do not mix polymer powder with another supplier monomer liquid and vice versa. Place the denture in water at 50 °C for 1 hour or at room temperature water for 24 hours prior placing in patient mouth.

Danger: Highly flammable liquid and vapour. / May cause respiratory irritation. / Causes skin irritation. / May cause an allergic skin reaction. / Causes serious eye irritation. / Keep away from heat/sparks/open flame/hot surfaces. - No smoking. / Keep container tightly closed. / Avoid breathing vapours. / Keep away from heat/sparks/open flames/hot surfaces. - No smoking / Keep container tightly closed. Avoid breathing vapours. / IF ON SKIN: Wash with plenty of soap and water. / Wear protective gloves/protective clothing/eye protection/face protection. / IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing./ If eye irritation persists: Get medical advice/attention. / Dispose of contents/containers in accordance with local regulation.

Personal protection: When working with product use protective nitrile gloves, when it is possible to inhale monomer use protective mask with filter A1, when processed polymerized product use protective mask with filter FFP1.

Side effects: Possibility of an allergic reaction to MMA and PMMA.

Storage:

- at room temperature
- in the place with good ventilation
- keep away from sunlight
- keep away from children

Notice: Any serious incident, that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

Warranty: Any technical consulting regarding our products, whether verbal, in writing or during demonstrations, is carried out to the best of our knowledge and is to be considered nonbinding. We have no control over use, application and processing, for which the user is responsible.

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