CC DISK PEEK

(SL) Navodila za uporabo / (EN) Directions for use / (DE) Gebrauchsanweisung / (CS) Návod k použití
(HR) Uputa za uporabu / (ES) Instrucciones de uso / (RU) Инструкции по использованию

CS

INTERDENT s.r.o.
Foerstrova 12, Strašnice
CZ-10000 Praha
T: +420/274 783 114
F: +420/274 820 130
E: interdent@interdent.cz

SR

INTERDENT d.o.o.
Egipatska 22, lok 3
RS-11070 Novi Beograd
T/F: +381/11 217 53 74
E: office@interdent-bg.com

HR

INTERDENT d.o.o.
Vinogradski odvojak 2d
HR-10431 Sveta Nedelja
T: +385/1 3873 644
F: +385/1 38736 17
E: interdent@interdent.hr

SK

Interdent SK s.r.o.
Za dráhou 21
SK-902 01 Pezinok
T: 0903 418 001
E: interdent@interdent.sk

INTERDENT®

Proizvajalec / Producer / Hersteller
Interdent d.o.o. · SI - 3000 CELJE
Opekarniška cesta 26 · T: +386 (0)3 425-62-00
E: info@interdent.cc · www.interdent.cc

CE 0197

Verzija: 04/2024
Datum: 29.05.2024

Made in Slovenia

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Intended use: fixed and removable dental prosthesis, abutments (where virtual library is available)

Indication:

- 4-part bridge up to 2-pontics for posterior area
- 6-part up to 2-pontics bridge anterior area
- partially removable
- telescopic work
- individual temporary abutment and healing abutment in anterior region
- fixed prosthesis all on 4 over implants
- reinforcement of total dentures (in region of alveolar ridge)
- Maryland bridges, inlay supported pontic
- endo-crowns
- root posts

Contraindications:

- where minimal dimensions for prosthetic appliances cannot be achieved
- possible allergic reactions to product content
- bruxism

Intended users: Product is intended to be used by professionals – dental technicians.

Patient selection criteria: Totally or partially edentulous patients.

Side effects: Possible allergic reaction to PEEK material.

Requirements for dimensions:

Removable partial dentures: design of frame 1.5 times more massive than in alloy. 0.5 mm of undercut of clasp design. Clasp activation is not possible. Material is recommended in case of known allergies to alloys and as a temporary solution.

Fixed dental prosthesis:

- Crowns & Bridges
Thickness: min 0.6 mm, possible thinner part in cervical area to 0.3 mm
Basal part of crown must be out of PEEK not veneering material
Connector anterior: 11 - 12 mm²
Connector posterior: 14 mm²
Cut back technique and garlands are recommended

- Design of the frame must follow end shape of the teeth
- All on 4 rehabilitation on implants:
 - Occlusal-cervical height: min 5 mm
 - Buccolingual width anterior: min 4 mm
 - Buccolingual width in area of titanium sleeve abutment: 6 mm

Procedure:

Design preparation: Make the design on the base of the model scan or digital impression of statues in the patient's mouth. On the base of the design dimension, place the right thickness of the disc CAD/CAM milling machine. The design is transferred to CAM SW on the milling machine and prosthetic appliance is milled out.

Processing and polishing: Use separating discs to cut the sprues. Use cross-sectioned carbide burs for removing any excess material. Parts that will not be covered by composite should be polished to high gloss. Use silicone eraser for ceramic and goat brushes with using diamond polishing paste to polish to high gloss. You can polish larger objects with pumice and cotton-based brushes. Process and polish with minimal pressure and slow speed (max. 8000 rpm). Polish delicate and fragile parts carefully, with lower speed and with smaller polishing brushes. The use of a hand piece is recommended.

Veneering: Sand-blast the surface that you want to veneer with $110 \mu\text{m Al}_2\text{O}_3$ under a pressure of 2 - 3 bar, and clean with compressed air. Do not touch sand-blasted surfaces anymore and do not use steam pressure. Apply bond to the surface (like Visio.link, Bredent GmbH) according to the instructions for use. Use composite that is appropriate for PEEK, and work according to the instructions for use.

Safety precaution: Use a suction unit, wear gloves, goggles, and a protective mask with filter FFP1 when processing the material.

Storage: room temperature, dry, protected against direct sunlight!

Notice: Any serious incident, that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

Warranty! Whether given verbally, in writing, or by practical instructions, our recommendations for use are based upon our own experience and trials, and can only be considered as standard values. Our products are subject to further development. Therefore, alteration in the construction and composition of working instructions is reserved.

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