I-BOND EASY

 (SL) Navodila za uporabo /	 (EN) Directions for use	
 (DE) Gebrauchsanweisung /	 (ES) Instrucciones de uso /	 (CS) Návod k použití
 (HR) Uputa za uporabu /	 (RU) Инструкции по использованию	

<p>CS INTERDENT s.r.o. Foerstrova 12, Strašnice CZ-10000 Praga T: +420/274 783 114 F: +420/274 820 130 E:interdent@interdent.cz</p>	<p>SR INTERDENT d.o.o. Zemunska 22, lok 3 RS-11070 Novi Beograd T/F: +381/11 217 53 74 www.interdent-bg.com</p>	<p>SK INTERDENT SK s.r.o. Za dráhou 21 SK-902 01 Pezínok T: 0903 418 001 E: interdent@interdent.sk</p>
<p>HR INTERDENT d.o.o Vinogradski odvojak 2d HR-10431 Sveta Nedelja T: +385/1 3873 644 F: +385/1 38736 17 E: interdent@interdent.hr</p>	<p>DE Interdent d.o.o. · SI - 3000 CELJE Opekarniška cesta 26 · T: +386 (0)3 425-62-00 E: info@interdent.cc · www.interdent.cc</p>	<p>CE 0197</p> <p>Verzija: 05/2022 Datum: 15.12.2022 Made in Slovenia</p>
<p>INTERDENT[®]</p>		
<p>Proizvajalec / Producer / Hersteller Interdent d.o.o. · SI - 3000 CELJE Opekarniška cesta 26 · T: +386 (0)3 425-62-00 E: info@interdent.cc · www.interdent.cc</p>		

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Melting and Casting: For investing, use phosphate bonded investments for crown and bridge work (Intervest K+B Speed or Interfine K+B Speed). Preheat the investment to about 850 to 900°C. Hold the end temperature of investment for a minimum 30 minutes. Refer to manufacturer’s instructions for use of the casting machines. For I-BOND EASY use an individual ceramic crucible to prevent contamination with other alloys. Clean crucible after every use. When melting by induction heating, start casting as soon as the ingots have collapsed and oxide net cracks. For melting by flame, rotate the reductive zone of the flame around the ingots. Start casting as soon as the bath begins to vibrate. Allow the cylinder to cool down slowly to room temperature and then deflask.

Processing: Process the frame in one direction using a hard-metal milling tool. Remove any sharp edges.

Ceramic firing: Use ceramic, suitable for porcelain-to-metal restoration.

Please be careful about the cooling procedure. Cobalt-based alloys must be cooled down slowly to 600 °C after each firing step (dentin, correction, glaze).

- Sandblast the surface by use of a pencil-blaster. Blast all of the surfaces on which the porcelain is to be applied with Al₂O₃ – Interlox 250 µm and clean the surface (ultrasonically cleaning, under vapor (steam cleaner) or boiling water (diethyl acetate).
- If you decide to perform oxide firing to check the surface, perform it at 960°C in a vacuum for 5–10 min. After oxide firing, removal of oxides is required – use aluminium oxide 250 microns and degrease the surface again.
- For a good bond between metal and ceramic it is mandatory to use a “metal-ceramic bond” such as VITA NP BOND paste. It should be applied in a thin layer, with brush (as Wash opaque). For firing in the ceramic furnace, check the below programme.**
- When the “metal-ceramic bond” phase is finished, proceed with opaque layering, to cover the metal shade of the construction.

Firing program using Vita NP Bond:						
Vt. in	→	↗	↘	Temp.	→	VAC
°C	min	min	°C/min	cca °C	min	min
600	6.00	6.00	60	960	1.00	6.00

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Cobalt-based dental casting alloy intended for metal-ceramic dental restorations without beryllium, nickel, cadmium and lead, type 4.

I-BOND EASY is a biocompatible non-precious alloy. I-BOND EASY and fulfils the recommendations of the standards EN ISO 22674 and EN ISO 9693.

Composition	Mass [%]	Properties	
Co	62,5	Type	4
Cr	27,2	Density	8,2 g/cm ³
W	8,2	Solidus, liquidus temperature	1380 °C, 1450 °C
Si	1,7	Casting temperature	1480 °C
Mn	< 1,0	Vickers hardness	HV 10 249
		Coefficient of thermal expansion	25 - 500 °C 14,4 x 10 ⁻⁶ K ⁻¹
			20 - 600 °C 14,6 x 10 ⁻⁶ K ⁻¹
		Yield stress	Rp 0,2 380 MPa
		Elongation	A5 18,1 %
		Tensile strength	Rm 593 MPa
		E-modul	E 167 GPa

Indications: Used for multiple units fixed dental prostheses for ceramic firing, also for appliances with thin sections that are subject to very high forces, e.g. removable partial dentures, clasps, thin veneered single crowns, full arch fixed dental prostheses or those