LC INTERTRAY

REF 0714

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GB Direction for use

Draw with the pencil the shape of the tray on stone model. Use »Space holder« wax if needed. Insulate stone model with Isoacryl and let it dry. Put the LC Intertray on the model and shape it. Use the cuttings for the holder of the tray. The tray can be polymerized in units that radiate haloge light from 420 nm to 480 nm wavelength or UV light from 350 nm to 400 nm wavelength. Put the stone model with the tray in the unit and polymerize it for 2 minutes. Remove the tray from the model, turn it upside down and polymerize it again for 2 minutes. Longer polymerization time improve the strength of the material.

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GB Direction for use

Draw with the pencil the shape of the tray on stone model. Use »Space holder« wax if needed. Insulate stone model with Isoacryl and let it dry. Put the LC Intertray on the model and shape it. Use the cuttings for the holder of the tray. The tray can be polymerized in units that radiate halogen light from 420 nm to 480 nm wavelength or UV light from 350 nm to 400 nm wavelength. Put the stone model with the tray in the unit and polymerize it for 2 minutes. Remove the tray from the model, turn it upside down and polymerize it again for 2 minutes. Longer polymerization time improve the strength of the material.

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