SL

PALICE ZA LOTANJE

Material za lotanje dentalnih zlitin na osnovi kobalta in niklja.

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EN

BRAZING RODS

Brazing material for cobalt- and nickel-based dental alloys.

Contain no nickel, beryllium, cadmium and lead and fulfill the recommendations of the standard EN ISO 9333.

Mas composition in %: Co 63.3, Cr 28.2, Si 4.2, Mo 2.4, B 1.4, C
Working temperature: 1310 °C

Solidus, liquidus temperature: 1130 °C, 1250 °C

Indications: Suitable for brazing of partial dentures, crowns and bridges. Suitable for ceramic veneering.

Contraindications: in case of known allergies to alloy content

Intended users: Product is intended to be used by professionals – dental technicians.

Patient selection criteria: **Totally or partially edentulous patients.**

Brazing: The wideness of soldering gap max 0,2 mm. Grind the gap and clean it by sand – blasting with pure Al₂O₃ (REF 0405 Interlox 110). Fix parts with the fire resistant block for soldering. Apply high temperature flux (Interflux REF 0495) and pre - heat soldering block to 600 °C. Heat up the gap by flame and put inside soldering part (Intersolder) wet with Interflux. Heat up the surrounding of the soldering area by flame till the temperature that Intersolder flows into soldering gap.

Finishing: Sand-blast the visible surface of the frame with Interlox 50 µm and then polish it with suitable grinding and polishing instruments for nickel alloys. At the end use Universal polishing paste for Co-Cr-Mo alloys REF 460 and polish up to high gloss. Clean polished surface with vapour or boiling water.

Safety warnings: Use suction unit, wear gloves, goggles and protective mask with filter FFP2 when processing the material, to remove and protect yourself against dust. Reciprocal Actions! In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

Side effects: Consider allergic hypersensitivities to contents of the alloy.

Notice: Any serious incident, that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

Warranty! Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our product are subject to further development. Therefore alterations in construction and composition are reserved.

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