CC DISK Ti5

Potrebno upoštevati navodila za uporabo / Consult instructions for use / Bitte beachten Sie die Bedienungsanleitung / Potrebno poštivati uputstva za uporabu / Viz návod k použití/ Consultar las instrucciones de uso / Dodrżujcie návod na použitie

Uporabno do / Use by / Gültig bis/ Upotrijebiti do/ Použítí / Usar antes / Použite do

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EN

Titanium-based cast discs, intended for milling in CAD/CAM machines to produce single crowns, larger bridges and implant superstructures, type 4.

CC DISK Ti5 are made of titanium, aluminium and vanadium and are used for appliances that have to be rigid and tough. CC DISK Ti5 are free of nickel, beryllium, cadmium and lead, and meet the requirements of EN ISO 22674 and EN ISO 9693.

Composition	m %	Technical data	
Ti	89,8	Type	4
Al	6,0	Density	4,43 g/cm ³
V	4,0	Vickers's hardness	HV10 353
Fe, O		Yield stress	R _{p0,2} 780 MPa
		Tensile strength	R _m 860 MPa
		Coefficient of thermal expansion	25- 500°C 9,8 x 10 ⁻⁶ K ⁻¹
		Elongation after fracture	A _{15mm} 10 %

Indications: Used for multiple-unit fixed dental prostheses for ceramics firing, also for appliances with thin sections that are subject to very high forces, e.g. removable partial dentures, clasps, thin veneered single crowns, full arch fixed dental prostheses or those with small cross-sections, bars, attachments, implant retained superstructures.

Contraindications: In case of known allergies to alloy content.

Intended users: The product is intended to be used by professionals – dental technicians.

Patient selection criteria: Totally or partially edentulous patients.

Milling: Use suitable milling drills for titanium. Smoothen the surface of the frame with special cross-cut drills for titanium or separating discs, always in the same direction to avoid blistering of the ceramic. Clean the surface with vapor or boiling distilled water.

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Ceramic firing and finishing:

- Sand blast the surface of the frame under a pressure of 2–3 bar with 110 µm Al₂O₃ (REF 405- Interlox 110 µm).
- Clean the surface with vapor or boiling distilled water.
- Use ceramics suitable for titanium that meet the demand of the standard EN ISO 9693 (e.g. VITA LUMEX AC). Firing has to be done according to the instructions of the manufacturer of the ceramics.
- Remove ceramics mechanically. Hydrofluoric acid (HF) corrodes the metal.
- On the cervical parts of the crowns, remove the oxides with glass pearls. Use silicone rubbers and polishing pastes for finishing.

Safety precautions: Metal powder is harmful to health. Use a suction unit and protective masks with an FFP1 filter.

Side effects: Consider allergic reactions to the alloy composition.

Caution: During dry milling of titanium, chips and swarf can ignite and cause fire. Wet milling is recommended.

Reciprocal actions: In case of occlusal or approximal contact of different alloys, electrochemical reactions may occur in very rare instances.

Notice: Any serious incident that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

Warranty: Whether given verbally, in writing, or by practical instructions, our recommendations for use are based upon our own experience and trials, and can only be considered as reference values. Our products are subject to further development. Therefore, alterations in the structure and composition are reserved.

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