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## GB

Nickel-based dental casting alloy intended for metal-ceramic dental restorations without beryllium, cadmium and lead, Type 4.

I-BOND LO is biocompatible non-precious nickel-based casting alloy for metal-ceramic dental restorations. I-BOND LO contains no beryllium, cadmium and lead and fulfills the recommendations of the standards EN ISO 22674 and EN ISO 9693-1. The low Vicker’s hardness of 250 (HV 10) leads to good and easy milling and grinding behaviour and gives perfect polished surfaces.

| Composition | (Mass-%): | Properties  |
|-------------|-----------|---|
| Ni          | 64,4      | Type 4  |
| Cr          | 22        | Density 8,3 g/cm <sup>3</sup>   |
| Mo          | 10        | Solidus, liquidus temperature 1200 <span> </span> °C, 1280 <span> </span> °C                        |
| Si          | 2,1       | Casting temperature 1340 <span> </span> °C  |
| Nb, Mn, B   | < 1       | Vickers hardness HV 10 250  |
|             |           | Coefficient of thermal expansion 25 - 500 <span> </span> °C 13.9 x 10 <sup>-6</sup> K <sup>-1</sup> |
|             |           | 20 - 600 <span> </span> °C 14.1 x 10 <sup>-6</sup> K <sup>-1</sup>                                  |
|             |           | Yield stress Rp 0,2 427 MPa   |
|             |           | Modulus of elasticity E aprox. 144.000 MPa  |
|             |           | Elongation A5 8,5 <span> </span> %  |

Recommendations for Use

Wax-up:

Do the wax-up as usual. The wall thickness in wax must not be **thinner than 0,35 mm**. Lead wax sprues indirectly. For wax sprues use round wax wires with **Ø 2,5 mm** for single crowns and **Ø 2,5 – 3,0 mm** for bridges. For greater prosthetic appliances with more than 4 teeth use a distribution bar with **Ø 4,0 – 4,5 mm**.

Investing and Casting:

For investing use phosphate bonded investments for crown and bridge work – (**Invest K+B Speed or Interfine K+B Speed**). Preheat the investment to about 850 to 900 °C. Hold end temperature of investment for a minimum of **30 minutes**. Refer to manufacturer’s instructions for use for the casting machines.

For I-BOND LO use an individual ceramic crucible to prevent contamination with other alloys. Clean crucible after every use. When melting by induction heating start casting as soon as the ingots have collapsed and oxide net cracks. For melting by flame heat rotate the reductive zone of the flame around ingots. Start casting as soon as the bath begins to vibrate. Allow the cylinder slow air cool down to the ambient temperature and deflask.

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**Ceramic:**

- Use ceramic, suitable for porcelain to metal restoration. Sand blast the surface by use of a pencil-blaster. Blast all of the surface on which the porcelain is to be applied with **Al<sub>2</sub>O<sub>3</sub> - Interlox 250 μm** and clean the surface (under vapour or boiling water, ethyl acetate).
- If **oxide firing** is carried out to check the surface, it need to be put on **960 °C** at vacuum for **5 - 10 minutes**. After cooling the frame has to show light yellow grey colour. Oxide layer should be removed with **Al<sub>2</sub>O<sub>3</sub> - Interlox 250 μm**. The frame is cleaned under vapour or boiling water
- The opaque is then applied on the surface by a first thin wash firing followed by a second covering opaque.
- Firing has to be done according to the ceramic manufacturer’s instructions.
- After each firing step (dentine bake, build-up and glazing) cool down in a normal manner**.
- Remove ceramic mechanically. Hydrofluoric acid (HF) corrodes the metal.

**Finishing:**

After firing of the metal-ceramic prosthetic work sand blast the visible surface of the frame with **Al<sub>2</sub>O<sub>3</sub> - Interlox 50 μm** and then polish it with suitable grinding and polishing instruments for nickel alloys. At the end use **Universal polishing paste for Co-Cr-Mo alloys REF 460** and polish up to high gloss. Clean polished surface with vapour or boiling water.

**Soldering and Welding:**

Soldering before firing of the frame can be carried out with suitable base-metal-solder and high temperature flux. For welding with laser use suitable base-metal welding wires. The requirements of EN ISO 9333 must be followed.

**Final Safety Notes!**

Metal dusts are harmful. Use dust extractor, wear gloves, googles and protective mask with filter FFP2 during handling with alloy.

**Side effects:**

Consider allergic hypersensitivities to contents of the alloy.

**Reciprocal Actions!**

In case of occlusal or approximal contact of different alloys electrochemical based reactions may very rarely occur.

**Warranty!**

Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our product are subject to further development. Therefore alteration in construction and composition is reserved.

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