# INTERACRYL PLAST

IA/2463

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## INTERACRYL PLAST

GB

(Acrylic material for custom impression trays)

Interacryl plast is a cadmium-free, self-curing material based on methyl methacrylate. It is supplied as a polymer and monomer and is intended for fabricating custom impression trays.

### Indications:

- Fabricating custom impression trays.

### Contraindications:

- not for removable dentures
- not for orthodontic

### Directions for use:

#### Mixing ratio at room temperature:

**Mixing ratio: 33,6 g powder / 8 ml liquid (4,2 g powder / 1 ml liquid)**

First pour the liquid into a silicone mixing vessel, then add the required amount of powder. Mix the two components, using a spatula, until the dough is no longer sticky and separates cleanly from the sides of the mixing vessel. This should take about 90 seconds at room temperature (22°C). Then the dough can be applied directly to the plaster model - which should have been coated with Isoacryl in advance. During the plastic phase, the material can be processed for approximately 6 minutes. At room temperature, the tray material will cure 12 minutes after mixing. To ensure that the tray retains its shape optimally, it should be post-cured by immersing it in water at a temperature of 60°C for a period of 4 minutes.

### Caution:

Do not mix polymer powder with another supplier monomer liquid and vice versa.

**Danger:** Highly flammable liquid and vapour. / May cause respiratory irritation. / Causes skin irritation. / May cause an allergic skin reaction. / Causes serious eye irritation. / Keep away from heat/sparks/open flame/hot surfaces. - No smoking. / Keep container tightly closed. / Avoid breathing vapours. / Keep away from heat/sparks/open flames/hot surfaces. - No smoking / Keep container tightly closed. Avoid breathing vapours. / IF ON SKIN: Wash with plenty of soap and water. / Wear protective gloves/protective clothing/eye protection/face protection. / IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing./ If eye irritation persists: Get medical advice/attention. / Dispose of contents/containers in accordance with local regulation.

### Personal protection:

We recommend wearing half facemask with filter A1.

### Storage:

- at room temperature
- in the place with good ventilation
- keep away from sunlight
- keep away from children

### Warranty:

Any technical consulting regarding our products, whether verbal, in writing or during demonstrations, is carried out to the best of our knowledge and is to be considered nonbinding. We have no control over use, application and processing, for which the user is responsible.

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