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GB

Titanium casted discs, intended for milling in CAD/CAM machines to produce crowns, shorter range bridges and implant superstructures, type 3.

CC DISK Ti2 are made of titanium class 2 and are free of nickel, beryllium, cadmium and lead. CC DISK Ti2 meet the requirements of standards EN ISO 22674 and EN ISO 9693-1.

Composition	(Mass- %):	Technical data
Ti	> 99 %	Alloy type according EN ISO 22674
Fe, C		Density 4,51 g/cm ³
		Vicker's hardness HV 10 145
		0,2 % - elongation limit Rp 0,2 275 MPa (N/mm ²)
		Tensile strenght 450 MPa (N/mm ²)
		Coefficient of thermal expansion 25 - 500 °C 9,6 x 10 ⁻⁶ K ⁻¹
		Ductile yield A5 20 %

Milling: Use suitable milling drills for titan. Smooth the surface of the frame with special cross-cut drills for titan or separating discs, always in same direction to avoid blistering of ceramic.

Ceramic firing and finishing:

- Sand blast the surface of the frame under pressure 2-3 bars with 110 µ Al₂O₃ (Interalox).
- Clean the surface with vapour or diethyleter. Do not use acids like HF.
- Use ceramic suitable for titanium that meet the demand of the standard EN ISO 9693-1 (E.G. VITA TITANKERAMIK). Firing has to be done according to the ceramic manufacturer's instructions.
- On the cervical parts of the crowns, remove the oxides with glass pearls. Use silicone rubbers and polishing pastes for finishing.

Safety precaution! Metal powder is health harmful. Use suction unit and protective masks with filter FFP1.

Side effects: Consider allergic reaction to the alloy composition.

CAUTION! During dry milling of titanium, chips and swarfs can ignite themselves and cause fire. Wet milling is recommended.

Warranty! Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our product are subject to further development. Therefore alterations in construction and composition are reserved.

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