CC DISK Ti5

Potrebno upoštevati navodila za uporabo / Consult instructions for use / Bitte beachten Sie die Bedienungsanleitung / Viz návod k použití /
Potrebno poštivati uputstva za upotrebu / Consultar las instrucciones de uso /
Необходимо учитывать инструкции по применению

Uporabno do / Use by / Gültig bis / Použití / Upotrijebiti do / Usar antes /
Использовать до

Kontrolna številka / Batch number / Seriennummer / Seriové číslo /
Serijski broj / Numero de serie / Серийный номер

Kataloška številka / Catalogue number / Katalognummer / Katalogové číslo /
Kataložski broj / Numero de catalogo / Каталогный номер

Proizvajalec / Manufacturer / Hersteller / Výrobce / Proizvodač / Fabricante /
Производитель

Datum proizvodnje / Date of manufacture / Herstellungsdatum /
Datum výroby / Datum proizvodnje / Fecha de fabricacion / Дата изготовления

SI Navodila za uporabo / GB Directions for use
D Gebrauchsanweisung / CZ Návod k použití
HR Uputa za uporabu

HR

INTERDENT d.o.o.
Vinogradski odvojak 2d
HR-10431 Sveta Nedelja
T: +385/1 3873 644
F: +385/1 38736 17
E: interdent@interdent.hr

SRB

INTERDENT d.o.o.
Zemunska 22, lok 3
RS-11070 Novi Beograd
T/F: +381/11 217 53 74
www.interdent-bg.com

CZ

INTERDENT s.r.o.
Foerstrova 12, Strašnice
CZ-10000 Praga
T: +420/274 783 114
F: +420/274 820 130
E: interdent@interdent.cz

SK

INTERDENT SK, s.r.o.
Liščie údolie 57
SK-84231 Bratislava
T: 00421 2 5440554,
F: 00421 2 5440555
E: interdent@interdent.sk

INTERDENT®

Proizvajalec/Producer/Hersteller
Interdent d.o.o. · SI - 3000 CELJE · Opekarniška cesta 26
T: +386 (0)3 425-62-00 · F: +386 (0)3 425-62-02
E: info@interdent.cc · www.interdent.cc

CE 0197

Verzija: 03/2018
Datum: 26.2.2018
Made in Slovenia

SI

Uliti diski na osnovi titana, namenjeni rezkanju v CAD/CAM strojih za izdelavo posameznih kron, večjih mostičkov in nadgradenj nad implantati, tip 5.

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GB

Titanium-based casted discs, intended for milling in CAD/CAM machines to produce single crowns, larger bridges and implant superstructures, type 5.

CC DISK Ti5 are made of titanium class 5 and are used for appliances that have to be rigid and tough. CC DISK Ti5 are free of nickel, beryllium, cadmium and lead and meet the requirements of EN ISO 22674 and EN ISO 9693-1.

Composition	(Mass- %):	Technical data
Ti	89,8	Alloy type according EN ISO 22674 5
Al	6	Density 4,43 g/cm ³
V	4	Vicker's hardness HV 10 353
Fe	0,2 % - elongation limit	Rp 0,2 780 MPa (N/mm ²)
	Tensile strenght	860 MPa (N/mm ²)
	Coefficient of thermal expansion	25 - 500 °C 9,8 x 10 ⁻⁶ K ⁻¹
	Ductile yield	A5 10 %

Milling: Use suitable milling drills for titan. Smooth the surface of the frame with special cross-cut drills for titan or separating discs, always in same direction to avoid blistering of ceramic.

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Ceramic firing and finishing:

1. Sand blast the surface of the frame under pressure 2-3 bars with 110 μ Al₂O₃ (Interalox).
2. Clean the surface with vapour or diethyleter. Do not use acids like HF.
3. Use ceramic suitable for titanium that meet the demand of the standard EN ISO 9693-1 (E.G. VITA TITANKERAMIK). Firing has to be done according to the ceramic manufacturer's instructions.
4. On the cervical parts of the crowns, remove the oxides with glass pearls. Use silicone rubbers and polishing pastes for finishing.

Safety precaution! Metal powder is health harmful. Use suction unit and protective msk with filter FFP1.

Side effects: Consider allergic reaction to the alloy composition.

CAUTION! During dry milling of titanium, chips and swarfs can ignite themselves and cause fire. Wet milling is recommended.

Warranty! Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our product are subject to further development. Therefore alterations in construction and composition are reserved.

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