# INTERACRYL HOT

(SL) Navodila za uporabo / (EN) Directions for use

(DE) Gebrauchsanweisung / (HR) Uputa za korištenje

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**INTERACRYL HOT** (heat-polymerizable resin)

Interacryl hot is a cadmium-free, heat-polymerizable methyl methacrylate-based resin, supplied in the form of a polymer (powder) and monomer (liquid). It is in accordance with demands of EN ISO 20795 - 1.

**Indications:**

- full dentures
- partial dentures

**Contraindications:** The possibility of an allergic reaction to methyl PMMA/MMA. Unpolymerised or poorly polymerised material don't place in the mouth

**Intended users:** Product is intended to be used by professionals – dental technicians.

**Patient selection criteria:** Totally or partially edentulous patients.

**Directions for use: (Preparation)** Invest the wax denture in the stone plaster. Carefully remove all the wax and treat the plaster model with an alginate liquid Isoacryl.

**Mixing ratio at room temperature:**

**Colours 3, 5, 10, 15, transparent: 23 g powder / 10 ml liquid (2,3 g powder / 1 ml liquid)**

Pour liquid in mixing bowl and then add powder. Start mixing immediately. Mix together thoroughly for 30 seconds. Leave the mixture to set for 10 minutes. When the mixture no longer sticks to sides of the mixing vessel, begin packing the mouldings. At an ambient temperature of approx. 20 °C the mixture will remain malleable for approx. 15 minutes. Pack the mouldings with a slight excess of mixture and gradually press them together pausing in between. The flasks should remain at room temperature at this time. To start the curing process, immerse the camped flasks in boiling water and leave to boil for 20 minutes. To cure extra thick dentures, after immersing in boiling water, turn of the heat and leave it to stand for 15 minutes. Then reheat and boil for 20 minutes. We recomend that the dentures be allowed to cool down as slowly as possible. Quenching the dentures in cold water can result in excessive shrinkage.

**Doughing time at 20 °C:**

- Colour 3, 5, 15: 13 ± 2 min (3, 5, 15)
- Colour 10, transparent: 7 ± 2 min

**Working time at 22 °C:** 15 minutes

**Curing time:** 20 minutes

**Caution:** Do not mix polymer powder with another supplier monomer liquid and vice versa. Place the appliance in water at 50 °C for 1 hour or at room temperature water for 24 hours prior placing in patient mouth.

**Danger:** Highly flammable liquid and vapour. / May cause respiratory irritation. / Causes skin irritation. / May cause an allergic skin reaction. / Causes serious eye irritation. / Keep away from heat/sparks/open flame/hot surfaces. - No smoking. / Keep container tightly closed. / Avoid breathing vapours. / Keep away from heat/sparks/open flames/hot surfaces. - No smoking / Keep container tightly closed. Avoid breathing vapours. / IF ON SKIN: Wash with plenty of soap and water. / Wear protective gloves/protective clothing/eye protection/face protection. / IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing./ If eye irritation persists: Get medical advice/attention. / Dispose of contents/containers in accordance with local regulation.

**Personal protection:** When working with product use protective nitrile gloves, when it is possible to inhale monomer use protective mask with filter A1, when processed polimerized product use protective mask with filter FFP1.

**Side effects:** Possibility of an allergic reaction to MMA and PMMA.

**Storage:**

- at room temperature
- in the place with good ventilation
- keep away from sunlight
- keep away from children

**Notice:** Any serious incident, that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

**Warranty:** Any technical consulting regarding our products, wheather verbal, in writing or during demonstrations, is carried out to the best of our knowledge and is to be considered nonbinding. We have no control over use, application and processing, for which the user is responsible.

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