Potrebno upoštevati navodila za uporabo / Consult instructions for use / Bitte beachten Sie die Bedienungsanleitung / Viz návod k použití / Potrebno poštivati uputstva za upotrebu / Consultar las instrucciones de uso / Необходимо учитывать инструкции по применению

Uporabno do / Use by / Gültig bis / Použiti / Upotrijebiti do / Usar antes / Использовать до

Kontrolna številka / Batch number / Seriennummer / Seriové číslo / Serijski broj / Numero de serie / Серийный номер

Kataloška številka / Catalogue number / Katalognummer / Katalogové číslo / Kataloški broj / Numero de catalogo / Каталожный номер

Medicinski pripomoček / Medical device / Medizinprodukt / Medicinski proizvod / Zdravotnički prostriedek / Producto sanitaria / Медицинское изделие

Edinstveni identifikator pripomočka / Unique device identifier / Einmalige Produktkennung / Jedinstvena identifikacija proizvoda / Уникальный идентификатор изделия / Jedinečným identifikátorem prostředku / Identificador único del producto

Proizvajalec / Manufacturer / Hersteller / Výrobce / Proizvodač / Fabricante / Производитель

Datum proizvodnje / Date of manufacture / Herstellungsdatum / Datum výroby / Datum proizvodnje / Fecha de fabricacion / Дата изготовления

Navodila za uporabo /

Directions for use

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Uputa za uporabu /

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EN

Cobalt-based casted discs, intended for producing prosthetic substitutes in CAD/CAM milling machines for metal-ceramic dental restorations, type 4.

CC DISK EASY CoCr is nickel, beryllium, cadmium and lead free and fulfills the recommendation of the standard EN ISO 22674 and EN ISO 9693 for dental alloys.

The alloy has excellent corrosion resistance and good milling properties. It is easy to drill in and gives perfectly polished surface. Only few oxides are formed during ceramic firing.

Composition	(%):	Technical data
Co	62,5	Type 4
Cr	27,2	Density 8,6 g/cm ³
W	8,2	Coefficient of thermal expansion 25 - 500 °C 15,4 x 10 ⁻⁶ K ⁻¹
Si	1,7	20 - 600 °C 15,9 x 10 ⁻⁶ K ⁻¹
Mn	<1,0	Elongation limit Rp 0,2 380 MPa
		Tensile strength Rm 539 MPa
		Modul of elasticity (GPa) 167 GPa
		Ductile yield A5 15,2 %

Indications: Used for multiple units fixed dental prostheses for ceramic firing, also for – for appliances with thin sections that are subject to very high forces, e.g. removable partial dentures, clasps, thin veneered single crowns, full arch fixed dental prostheses or those with small cross-sections, bars, attachments, implant retained superstructures.

Contraindications: in case of known allergies to alloy content

Intended users: Product is intended to be used by professionals – dental technicians.

Patient selection criteria: Totally or partially edentulous patients.

Ceramic firing: You can use all types of ceramic which you also use in conventional method on cast crowns and bridges (VITA, Ivoclar, Ducera, etc.).

Please be careful about cooling procedure. Cobalt-based alloys must be cooled down slowly to 600 °C after each firing step (dentin, correction, glaze).

- Remove connectors and smoothen the sharp edges. Sand-blast the surface with 250 microns aluminum oxide on 2-3 bar air pressure (REF 404 Interlox 250 micr).

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