Potrebno upoštevati navodila za uporabo / Consult instructions for use / Bitte beachten Sie die Bedienungsanleitung / Viz návod k použití / Potrebno poštivati uputstva za upotrebu / Consultar las instrucciones de uso / Необходимо учитывать инструкции по применению

Uporabno do / Use by / Gültig bis / Použiti / Upotrijebiti do / Usar antes / Использовать до

Kontrolna številka / Batch number / Seriennummer / Seriové číslo / Serijski broj / Numero de serie / Серийный номер

Katalogska številka / Catalogue number / Katalognummer / Katalogové číslo / Kataloški broj / Numero de catalogo / Каталожный номер

Medicinski pripomoček / Medical device / Medizinprodukt / Medicinski proizvod / Zdravotnički prostriedek / Producto sanitaria / Медицинское изделие

Edinstveni identifikator pripomočka / Unique device identifier / Einmalige Produktkennung / Jedinstvena identifikacija proizvoda / Уникальный идентификатор изделия / Jedinečným identifikátorem prostředku / Identificador único del producto

Proizvajalec / Manufacturer / Hersteller / Výrobce / Proizvodač / Fabricante / Производитель

Datum proizvodnje / Date of manufacture / Herstellungsdatum / Datum výroby / Datum proizvodnje / Fecha de fabricacion / Дата изготовления

INTERDENT®

Proizvajalec / Producer / Hersteller
Interdent d.o.o. · SI - 3000 CELJE
Opekarniška cesta 26 · T: +386 (0)3 425-62-00
E: info@interdent.cc · www.interdent.cc

Navodila za uporabo /

Directions for use

Gebruiksaanweisung /

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Uputa za uporabu /

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EN

Cobalt-based casted discs, intended for producing prosthetic substitutes in CAD/CAM milling machines for metal-ceramic dental restorations, type 4.

CC DISK EASY CoCr is nickel, beryllium, cadmium and lead free and fulfills the recommendation of the standard EN ISO 22674 and EN ISO 9693 for dental alloys.

The alloy has excellent corrosion resistance and good milling properties. It is easy to drill in and gives perfectly polished surface. Only few oxides are formed during ceramic firing.

Composition	(%):	Technical data
Co	62,5	Type
Cr	27,2	Density
W	8,2	Temperature solidus, liquidus
Si	1,7	Vickers hardness
Mn	<1,0	Coefficient of thermal expansion
		HV 10
		249
		25 - 500 °C
		14,4 x 10 ⁻⁶ K ⁻¹
		20 - 600 °C
		14,6 x 10 ⁻⁶ K ⁻¹
		Elongation limit
		Rp 0,2
		380 MPa
		Tensile strength
		Rm
		539 MPa
		Modul of elasticity (GPa)
		167 GPa
		Ductile yield
		A5
		15,2 %

Indications: Used for multiple units fixed dental prostheses for ceramic firing, also for – for appliances with thin sections that are subject to very high forces, e.g. removable partial dentures, clasps, thin veneered single crowns, full arch fixed dental prostheses or those with small cross-sections, bars, attachments, implant retained superstructures.

Contraindications: in case of known allergies to alloy content

Intended users: Product is intended to be used by professionals – dental technicians.

Patient selection criteria: Totally or partially edentulous patients.

Ceramic firing: You can use all types of ceramic which you also use in conventional method on cast crowns and bridges (VITA, Ivoclar, Ducera, etc.).

Please be careful about cooling procedure. Cobalt-based alloys must be cooled down slowly to 600 °C after each firing step (dentin, correction, glaze).

- Remove connectors and smoothen the sharp edges. Sand-blast the surface with 250 microns aluminum oxide on 2-3 bar air pressure (REF 404 Interlox 250 micr).

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-
- Clean the surface in distilled water with ultrasound, degrease it with ethyl acetate or under vapor.
- Carry out oxide firing (only to control the surface) 5 min at 960°C. After oxide firing removing of oxides is required - use aluminum oxide 250 microns and degrease the surface again.
- For a good bond between metal and ceramic it is mandatory to use a "metal-ceramic bond". VITA NP BOND paste should be applied in a thin layer, with brush (as Wash opaque). For firing in the ceramic furnace check the below program.**

Firing program using Vita NP Bond:						
Vt. in	→	↗	↗	Temp.	→	VAC
°C	min	min	°C/min	cca °C	min	min
600	6.00	6.00	60	960	1.00	6.00

-
-
-
-
- When "metal-ceramic bond" phase is finished, proceed with opaque layering, to cover the metal shade of the construction.
- For firing procedure, consult ceramic manufacturer's instructions for use.**
- After each firing step (dentine bake, correction bake, glaze firing), cool down slowly to 600 °C.

Finishing: After last firing step finish the visible surface of the alloy with polishing rubbers and polishing paste (REF 461 Diamond polishing paste) to high gloss. On the cervical parts of the crowns, remove the oxides with glass pearls or very carefully with 50 microns aluminum oxide (REF 406 Interlox 50 micr.) Be careful not to damage the gingival fit!

Safety warnings: Use suction unit, wear gloves, googles and protective mask with filter FFP2 when processing the material, to remove and protect yourself against dust.

Reciprocal Actions! In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

Side effects: Consider allergic hypersensitivities to contents of the alloy.

Notice: Any serious incident, that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

Warranty: Technical data based on our own researches and knowledge and can be used as orientation values. The user itself is responsible for appropriate use of the product.

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