CC DISK NF CoCr

(SL) Navodila za uporabo / (EN) Directions for use / (DE) Gebrauchsanweisung
(HR) Uputa za uporabu / (CS) Návod k použití / (RU) Инструкции по использованию /
(KK) Пайдалану бойынша нұсқаулық / (ES) Instrucciones de uso / (FR) Instructions de procédure
(SK) Návod na použitie / (UK) Інструкція для застосування

CS

INTERDENT s.r.o.
Foerstrova 12, Strašnice
CZ-10000 Praga
T: +420/274 783 114
F: +420/274 820 130
E: interdent@interdent.cz

SK

Interdent SK s.r.o.
Za dráhou 21
SK-902 01 Pezinok
T: 0903 418 001
E: interdent@interdent.sk

HR

INTERDENT d.o.o
Vinogradski odvojak 2d
HR-10431 Sveta Nedelja
T: +385/1 3873 644
F: +385/1 38736 17
E: interdent@interdent.hr

SR

INTERDENT d.o.o.
Egipatska 22, lok 3
RS-11070 Novi Beograd
T/F: +381/11 217 53 74
www.interdent-bg.com

Interdent d.o.o. · SI - 3000 CELJE
Opekarniška cesta 26 · T: +386 (0)3 425-62-00
E: info@interdent.cc · www.interdent.cc

CE 0197

Version: 13/2026
Date: 20.2.2026
Made in Slovenia

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Cobalt-based cast discs intended for producing prosthetic substitutes in CAD/CAM milling machines for metal-ceramic dental restorations, type 4.

CC DISK NF CoCr is nickel-, beryllium-, cadmium- and lead-free, and meets the recommendations of standards EN ISO 22674 and EN ISO 9693 for dental alloys. The Vickers hardness 281 (HV10) value leads to good milling and drilling behaviour and gives a perfectly polished surface. Only few oxides are formed during ceramic firing.

Composition	m %	Technical data	
Co	63,0	Type	4
Cr	24,0	Density	8,3 g/cm ³
W	8,0	Coefficient of thermal expansion	25-500 °C
Mo	3,0		20-600 °C
Si	1,0	Vickers hardness	HV 10
Nb	<1,0	Yield stress	R _{p0,2}
		Tensile strength	R _m
		Modul of elasticity	E
		Elongation after fracture	A _{15mm}

Indications: Used for multiple-unit fixed dental prostheses for ceramics firing, also for appliances with thin sections that are subject to very high forces, e.g. removable partial dentures, clasps, thin veneered single crowns, full arch fixed dental prostheses or those with small cross-sections, bars, attachments, implant retained superstructures.

Contraindications: In case of known allergies to alloy content.

Intended users: The product is intended to be used by professionals– dental technicians.

Patient selection criteria: Totally or partially edentulous patients.

Ceramic firing: Use ceramics suitable for porcelain to metal restoration, with an appropriate coefficient of thermal expansion (VITA, Ivoclar, Ducera, etc.).

Please be careful about the cooling procedure. Cobalt-based alloys must be cooled down slowly to 600 °C after each firing step (dentin, glaze, correction).

1. Remove connectors and smoothen the sharp edges. For surface finishing, use carbide burs. Sand-blast the surface with 250-micron aluminium oxide at 2–3 bar air pressure (REF 404 Interlox 250 micr).
2. Clean the surface with steam or boiling distilled water. Do not touch the surface with your hands afterwards.
3. If oxide firing is carried out to check the surface, it needs to be fired at 960 °C in a vacuum for 5–10 minutes. The oxide layer should be removed with Al₂O₃ – REF 404 Interlox 250 µm at 2–3 bar of pressure. The frame is cleaned under vapour or boiling distilled water.
4. For the firing procedure, consult the ceramic manufacturer's instructions for use.
5. After each firing step (dentine bake, glaze firing, correction bake), cool down slowly to 600 °C.
6. Remove ceramics mechanically. Hydrofluoric acid (HF) corrodes the metal.

Finishing: After the last firing step, finish the visible surface of the alloy with polishing rubbers and a polishing paste for Co-Cr alloys (REF 461- Diamond polishing paste) to high gloss. Very carefully remove oxides with glass pearls or with 50 µm aluminium oxide (REF 406- Interlox 50 µm) on the cervical parts of the crowns (ensure that gingival fit is not damaged).

Safety warnings: Use a suction unit, wear gloves, goggles and a protective mask with an FFP2 filter when processing the material to remove and protect yourself against dust.

Reciprocal actions: In case of occlusal or approximal contact of different alloys, electrochemical reactions may occur in very rare instances.

Side effects: Consider allergic hypersensitivities to the contents of the alloy.

Notice: Any serious incident that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

Warranty: Whether given verbally, in writing, or by practical instructions, our recommendations for use are based upon our own experience and trials, and can only be considered as reference values. Our products are subject to further development. Therefore, alterations in the structure and composition are reserved.

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