**CC DISK Zr**  
**CC DISK Zr HT**

**CC DISK Zr HT Multilayer**

**CC DISK Zr Multicolour**

Ⓘ Navodila za uporabo / Ⓒ Directions for use / Ⓓ Gebrauchsanweisung / Ⓒ Návod k použití  
Ⓗ Uputa za uporabu / Ⓔ Instrucciones de uso / Ⓓ Инструкции по использованию

Proizvajalec / Producer / Hersteller  
Interdent d.o.o. · SI - 3000 CELJE  
Opekarniška cesta 26 · T: +386 (0)3 425-62-00  
E: info@interdent.cc · [www.interdent.cc](http://www.interdent.cc)

CE 0197

Verzija: 05/2019  
Datum: 06.03.2019

Made in Slovenia

**CZ**

INTERDENT s.r.o.  
Foerstrova 12, Strašnice  
CZ-10000 Praga  
T: +420/274 783 114  
F: +420/274 820 130  
E:interdent@interdent.cz

**SRB**

INTERDENT d.o.o.  
Zemunska 22, lok 3  
RS-11070 Novi  
Beograd  
T/F: +381/11 217 53 74  
www.interdent-bg.com

**HR**

INTERDENT d.o.o  
Vinogradski odvojak 2d  
HR-10431 Sveta Nedelja  
T: +385/1 3873 644  
F: +385/1 38736 17  
E: interdent@interdent.hr

**SK**

INTERDENT SK, s.r.o.  
Líščie údolie 57  
SK-84231 Bratislava  
T: 00421 2 5440554,  
F: 00421 2 5440555  
E: interdent@interdent.sk

SI Navodila za uporabo .....	4, 5, 6, 7
GB Directions for use .....	8, 9, 10, 11
D Gebrauchsanweisung .....	12, 13, 14, 15
CZ Návod k použití .....	16, 17, 18, 19
ES Instrucciones de uso .....	20, 21, 22, 23
HR Uputa za uporabu .....	24, 25, 26, 27
RUS Инструкции по использованию .....	28, 29, 30, 31

CC DISK Zr (beli in obarvani, opačni) CC DISK Zr HT (brezbarvni in obarvani, super translucetni), CC DISK Zr HT Multilayer (večslojni obarvani, super translucetni) in CC DISK Zr Multicolour (stopenjsko obarvani, translucetni) so diski iz biokompatibilnega predsintiranega cirkonijevega dioksida, stabiliziranega z itrijevim oksidom. Uporabljajo se v CAD/CAM rezkalnih strojih. CC DISK Zr (HT, HT Multilayer, Multicolour) ustrezajo zahtevam standarda EN ISO 6872, tip II, razred 5.

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**CC DISK Zr (white and coloured, opaque), CC DISK Zr HT (pure and coloured, high translucent), CC DISK Zr HT Multilayer (multilayer coloured, high translucent), CC DISK Zr Multicolour (gradually coloured, translucent)** are discs made out of biocompatible pre-sintered zirconium dioxide stabilized with yttrium oxide. They are used in in CAD/CAM milling machines. CC DISK Zr (HT, HT Multilayer, Multicolour) comply with demands of the standard EN ISO 6872, type II, class 5.

**Indications:** CC DISK Zr, CC DISK Zr HT, CC DISK Zr HT Multilayer and CC DISK Zr Multicolour are used for the production of full anatomical restorations, as for classical frames, meant for porcelain veneering (like caps, frames and all ceramic restorations with max. 2 pontics in anterior and 4 pontics in posterior section), telescopes, bar constructions, Maryland bridges, inlays, onlays and individual implant abutments.

**Contraindications:** Using in case of bruxism, inadequate strength of tooth structure or insufficient space is not recommended. Usage is also not recommended in case of inadequate oral hygiene or inadequate coordination of masticatory movements. Prosthetic substitutes made out of discs must not be used in case of unsatisfactory preparation of teeth and in case of mistakes in prosthetic substitute.

Physical and chemical characteristics	CC DISK Zr	CC DISK Zr coloured	CC DISK Zr HT	CC DISK Zr HT coloured	CC DISK Zr Multilayer	CC DISK Zr Multicolour coloured
Colour	White opaque	A1, A2, A3 opaque	White translucent	A1, A2, A3 high translucent	Multilayer coloured: A1, A2, A3, A3.5, B1, B2, B3, C1, C2, D2, D3, high translucent	Gradually coloured: A2 (A1 - A2,5), A3 (A2 - A3,5) translucent
Density [g/cm <sup>3</sup> ]	>6,05	>6,00	>6,00	>6,00	>6,02	> 6,00
Flexural strength [MPa]	1200±200	1000±200	1000±200	1000±200	≥ 800	1000±200
Thermal expansion [10 <sup>-6</sup> K <sup>-1</sup> ]	10,6	10,6	10,7	10,7	10,5	10,7
Translucency [%]	/	/	41	41	46	41
Radioactivity [Bq/g]	< 0,2	< 0,2	< 0,2	< 0,2	< 0,1	< 0,2

**Demands for preparation:** Doctor of dental medicine is responsible for proper tooth preparation. Please consider that the connector section has to be dimensioned larger, depending on the construction.

**Below written minimal wall thickness dimension are reference values:**

Posterior section	Occlusal 0,9 mm	Circular 0,5 mm
Anterior section	Incisal/occlusal 0,7mm	Circular 0,5 mm
Abutment crowns of bridge frameworks with 2 pontics	Incisal/occlusal 1,0mm	Circular 0,7 mm

**The connector surface should be at least:**

Anterior bridge	With one pontic	7 mm <sup>2</sup>
Anterior bridge	With two pontics	9 mm <sup>2</sup>
Posterior bridge	With one pontic	9 mm <sup>2</sup>
Posterior bridge	With two pontics	12 mm <sup>2</sup>

**Milling:** Milling can be done in any dental CAD/CAM machine. Milling machine's manufacturer must inform user about suitable cutters and revolutions per minute for milling CC DISK Zr. Discs must not come into contact with liquid or oil during milling in order to avoid strength loose or staining.

Wet milling of presintered discs is possible but construction must be dried out completely before sintering process. Drying can be done in drying chamber on 80 °C or under infrared light. As an alternative low-temperature firing program can also be selected, which does not harm the construction. Notches, sharp edges or slight damages, also in inner part, can reduce sintered construction strength, therefore it is important that rounded and smooth surface is ensured after milling.

**Colouring:** If required, the frame may be coloured to desired colour with standard coloring solutions that are suitable for zirconium oxide ceramics. In order to achieve uniform result, it must be ensured that the framework is completely clean, greaseless and dry before coloring. For the correct working procedure, instructions of the coloring solutions manufacturer should be followed.

## Sintering

## Recommended sintering program

## CC DISK Zr, CC DISK Zr HT, CC DISK Zr Multicolor

	Standard program	Fast program (single crowns/ bridges with three pontics)
Heating rate	6 °C per minute up to 700 °C	45 °C per minute up to 700 °C
Maintenance time	without	without
Heating rate	5 °C per minute up to 1450 °C	30 °C per minute up to 1450 °C
Final temperature	1450 °C	1450 °C
Maintenance time at the final temperature	120 minutes	50 minutes
Cooling	in a closed furnace	in a closed furnace

Fast program is intended to be used for single crowns and bridges with maximum 3 pontics.

In order to achieve the maximum predicted transparency of the material, it is recommended to use only a standard program.

## CC DISK Zr HT Multilayer

	1-7 units	Above 7 units
Heating rate	8 °C per minute up to 300 °C	6 °C per minute up to 300 °C
Maintenance time	without	without
Heating rate	5 °C per minute up to 1520 °C	5 °C per minute up to 1520 °C
Final temperature	1520 °C	1520 °C
Maintenance time at the final temperature	90 minutes	90 minutes
Cooling	controlled 4,5 °C / min to 800 °C, then cooling in a closed furnace	controlled 3 °C / min to 800 °C, then cooling in a closed furnace

Sintering can be done in all standard furnaces that are suitable for zircon oxide ceramic. Follow the instruction for use of the furnace.

Enlargement factor depends on the lot and it is written on the edge of the disc.

We recommend sintering construction in sintering crucible, placed in suitable granular material. One of the suitable granular type are ZrO<sub>2</sub> sintering spheres. If you use crucible with granular material, place the construction with the occlusal side down on the granular material in the crucible. Prevent contacts with other frames or sintering carriers. If you use sintering carriers or wires than make sure that they do not come into contact with the walls of the caps.

Dental constructions and crowns made out of zirconium oxide can be sintered only in suitable furnaces. If other ceramic materials are fired or alloy is sintered in the same furnace, this can influence the sintering behavior or can cause local discoloration.

**Ceramic veneering:** You can use any ceramic suitable for veneering zirconium oxide ceramic. Follow the instruction for use of ceramic for veneering.

**Finishing after final sintering:** Any finishing after final sintering must be done with water spray to prevent local overheating or damage the material.

**Storing:** If discs are stored in original packaging (cool and dry, meaning 10 °C - 35 °C and 30 % - 70 % relative humidity) than changes in characteristics of the material are not expected, also if stored for a longer period. Do not expose discs to direct sunlight.

**Disposal of the residual material:** Residual part of discs and sintered part of discs are chemically inert. Dispose them in accordance with local regulation.

**Safety precaution!** Powder is harmful to health. Use suction unit. At manual treatment use protective mask FFP1 and wear goggles.

**Side effects:** Consider possibility of allergic hypersensitivity to the powder.

**Warranty!** Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our product is subject to further development. Therefore, alteration in construction and composition of working instructions is reserved.

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